

LOT		MX 407 L		Rev 1		CONTROL		L		 (Exp.)		2017-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.9	2.0	2.1	2.1	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.38	2.41	2.40	2.40						± 0.15	
HB HGB	g/dl	6.2	6.2	6.2	6.1	6.1						± 0.4	
	g/l	62	62	62	61	61						± 4	
	mmol/l	3.85	3.85	3.85	3.79	3.79						± 0.25	
HT HCT	%	16.8	16.7	16.9	16.8	16.8						± 2.0	
	l/l	0.168	0.167	0.169	0.168	0.168						± 0.020	
VGM MCV	µm ³ ; fl	70	70	70	70	70						± 4	
TGMH MCH	pg	25.8	26.1	25.7	25.4	25.4						± 2.0	
	fmol	1.60	1.62	1.60	1.58	1.58						± 0.12	
CCMH MCHC	g/dl	36.9	37.2	36.8	36.3	36.3						± 3.0	
	g/l	369	372	368	363	363						± 30	
	mmol/l	22.92	23.11	22.82	22.55	22.55						± 1.86	
IDR RDW	%	14.0	14.0	N/A	14.5	13.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	75	72	68	74	75						± 20	
VPM MPV	µm ³ ; fl	9.2	9.0	N/A	9.0	8.3						± 2.0	
LYMPHO	#	1.24	1.32	N/A	1.16	1.28						± 0.40	
	%	65.0	66.0	N/A	55.0	64.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.29	0.20						± 0.20	
	%	9.0	10.0	N/A	14.0	10.0						± 6.0	
GRANULO	#	0.49	0.48	N/A	0.65	0.52						± 0.40	
	%	26.0	24.0	N/A	31.0	26.0						± 7.0	

LOT		MX 407 N		CONTROL N				 (Exp.)		2017-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.6	7.6	7.7	7.5	7.6					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.57	4.60	4.60	4.60	4.61					± 0.18
HB HGB	g/dl	13.4	13.3	13.3	13.3	13.3					± 0.5
	g/l	134	133	133	133	133					± 5
	mmol/l	8.32	8.26	8.26	8.26	8.26					± 0.31
HT HCT	%	36.6	36.8	37.3	36.8	36.9					± 2.5
	l/l	0.366	0.368	0.373	0.368	0.369					± 0.025
VGM MCV	µm ³ ; fl	80	80	81	80	80					± 4
TGMH MCH	pg	29.3	28.9	28.9	28.9	28.9					± 2.0
	fmol	1.82	1.80	1.80	1.80	1.79					± 0.12
CCMH MCHC	g/dl	36.7	36.1	35.7	36.1	36.1					± 3.0
	g/l	367	361	357	361	361					± 30
	mmol/l	22.76	22.44	22.17	22.44	22.40					± 1.86
IDR RDW	%	13.5	13.5	N/A	14.0	12.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	257	255	250	254	253					± 40
VPM MPV	µm ³ ; fl	8.5	8.3	N/A	8.6	7.7					± 2.0
LYMPHO	#	2.66	2.66	N/A	2.25	2.66					± 0.40
	%	35.0	35.0	N/A	30.0	35.0					± 5.0
MONO	#	0.46	0.53	N/A	0.53	0.49					± 0.30
	%	6.0	7.0	N/A	7.0	6.5					± 4.0
GRANULO	#	4.48	4.41	N/A	4.73	4.45					± 0.60
	%	59.0	58.0	N/A	63.0	58.5					± 6.0

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B



LOT		MX 407 H		CONTROL		H		 (Exp.)		2017-12-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE								
		MICROS		MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60														
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.4	20.1	20.3	20.1	20.3													± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.77	5.65	5.76	5.73	5.71													± 0.20	
HB HGB	g/dl	17.9	17.7	17.7	17.9	17.8													± 0.6	
	g/l	179	177	177	179	178													± 6	
	mmol/l	11.12	10.99	10.99	11.12	11.05													± 0.37	
HT HCT	%	50.8	49.7	50.7	50.4	50.2													± 3.0	
	l/l	0.508	0.497	0.507	0.504	0.502													± 0.030	
VGM MCV	µm ³ ; fl	88	88	88	88	88													± 4	
TGMH MCH	pg	31.0	31.3	30.7	31.2	31.2													± 2.0	
	fmol	1.93	1.95	1.91	1.94	1.94													± 0.12	
CCMH MCHC	g/dl	35.3	35.6	34.9	35.5	35.4													± 3.0	
	g/l	353	356	349	355	354													± 30	
	mmol/l	21.89	22.11	21.69	22.04	22.00													± 1.86	
IDR RDW	%	13.0	13.0	N/A	13.5	12.0													± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	525	515	525	530	510													± 45	
VPM MPV	µm ³ ; fl	8.3	8.2	N/A	8.5	7.6													± 2.0	
LYMPHO	#	3.26	3.02	N/A	2.71	3.15													± 1.00	
	%	16.0	15.0	N/A	13.5	15.5													± 5.0	
MONO	#	1.12	1.21	N/A	0.80	0.91													± 0.80	
	%	5.5	6.0	N/A	4.0	4.5													± 4.0	
GRANULO	#	16.01	15.88	N/A	16.58	16.24													± 1.40	
	%	78.5	79.0	N/A	82.5	80.0													± 7.0	