

| LOT                      |                                                        | MX 430 L        |  | CONTROL  |  | L      |  |  (Exp.) |  | 2021-10-05<br>( YYYY - MM - DD ) |  |                         |        |  |        |         |        |  |
|--------------------------|--------------------------------------------------------|-----------------|--|----------|--|--------|--|--------------------------------------------------------------------------------------------|--|----------------------------------|--|-------------------------|--------|--|--------|---------|--------|--|
|                          |                                                        | Rev 1           |  | ABX Lyse |  |        |  |                                                                                            |  |                                  |  |                         |        |  |        |         |        |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        |                 |  |          |  |        |  |                                                                                            |  |                                  |  | TOLERANCES<br>TOLERANCE |        |  |        |         |        |  |
|                          |                                                        | MICROS          |  | MICROS   |  | MICROS |  | MICROS                                                                                     |  | ARGOS - HELIOS                   |  |                         | PENTRA |  | PENTRA |         | PENTRA |  |
|                          |                                                        | ES60<br>Care ST |  | CRP 200  |  |        |  | 60 8P - 60 18P<br>ADVIA 60                                                                 |  |                                  |  |                         |        |  |        |         |        |  |
| GB WBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 2.0             |  | 2.0      |  |        |  | 2.0                                                                                        |  |                                  |  |                         |        |  |        | ± 0.4   |        |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 2.35            |  | 2.30     |  |        |  | 2.35                                                                                       |  |                                  |  |                         |        |  |        | ± 0.15  |        |  |
| HB HGB                   | g/dl                                                   | 6.0             |  | 6.0      |  |        |  | 6.0                                                                                        |  |                                  |  |                         |        |  |        | ± 0.4   |        |  |
|                          | g/l                                                    | 60              |  | 60       |  |        |  | 60                                                                                         |  |                                  |  |                         |        |  |        | ± 4     |        |  |
|                          | mmol/l                                                 | 3.73            |  | 3.73     |  |        |  | 3.73                                                                                       |  |                                  |  |                         |        |  |        | ± 0.25  |        |  |
| HT HCT                   | %                                                      | 16.5            |  | 16.1     |  |        |  | 16.5                                                                                       |  |                                  |  |                         |        |  |        | ± 2.0   |        |  |
|                          | l/l                                                    | 0.165           |  | 0.161    |  |        |  | 0.165                                                                                      |  |                                  |  |                         |        |  |        | ± 0.020 |        |  |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 70              |  | 70       |  |        |  | 70                                                                                         |  |                                  |  |                         |        |  |        | ± 4     |        |  |
| TGMH MCH                 | pg                                                     | 25.5            |  | 26.1     |  |        |  | 25.5                                                                                       |  |                                  |  |                         |        |  |        | ± 2.0   |        |  |
|                          | fmol                                                   | 1.59            |  | 1.62     |  |        |  | 1.59                                                                                       |  |                                  |  |                         |        |  |        | ± 0.12  |        |  |
| CCMH MCHC                | g/dl                                                   | 36.5            |  | 37.3     |  |        |  | 36.5                                                                                       |  |                                  |  |                         |        |  |        | ± 3.0   |        |  |
|                          | g/l                                                    | 365             |  | 373      |  |        |  | 365                                                                                        |  |                                  |  |                         |        |  |        | ± 30    |        |  |
|                          | mmol/l                                                 | 22.65           |  | 23.14    |  |        |  | 22.65                                                                                      |  |                                  |  |                         |        |  |        | ± 1.86  |        |  |
| IDR RDW                  | %                                                      | 14.5            |  | 14.0     |  |        |  | 13.5                                                                                       |  |                                  |  |                         |        |  |        | ± 3.0   |        |  |
| PLAQ. PLTS               | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 77              |  | 75       |  |        |  | 78                                                                                         |  |                                  |  |                         |        |  |        | ± 20    |        |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 9.3             |  | 9.3      |  |        |  | 8.8                                                                                        |  |                                  |  |                         |        |  |        | ± 2.0   |        |  |
| LYMPHO                   | #                                                      | 1.36            |  | 1.41     |  |        |  | 1.35                                                                                       |  |                                  |  |                         |        |  |        | ± 0.40  |        |  |
|                          | %                                                      | 68.0            |  | 70.5     |  |        |  | 67.5                                                                                       |  |                                  |  |                         |        |  |        | ± 8.0   |        |  |
| MONO                     | #                                                      | 0.20            |  | 0.20     |  |        |  | 0.20                                                                                       |  |                                  |  |                         |        |  |        | ± 0.20  |        |  |
|                          | %                                                      | 7.0             |  | 7.5      |  |        |  | 7.5                                                                                        |  |                                  |  |                         |        |  |        | ± 6.0   |        |  |
| GRANULO                  | #                                                      | 0.50            |  | 0.44     |  |        |  | 0.50                                                                                       |  |                                  |  |                         |        |  |        | ± 0.40  |        |  |
|                          | %                                                      | 25.0            |  | 22.0     |  |        |  | 25.0                                                                                       |  |                                  |  |                         |        |  |        | ± 7.0   |        |  |

| LOT                      |                                                        | MX 430 N        |         | Rev 1  |        | CONTROL N                  |                |  (Exp.) |        | 2021-10-05<br>( YYYY - MM - DD ) |                         |
|--------------------------|--------------------------------------------------------|-----------------|---------|--------|--------|----------------------------|----------------|--------------------------------------------------------------------------------------------|--------|----------------------------------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        | ABX Lyse        |         |        |        |                            |                |                                                                                            |        |                                  | TOLERANCES<br>TOLERANCE |
|                          |                                                        | MICROS          | MICROS  | MICROS | MICROS | MICROS                     | ARGOS - HELIOS | PENTRA                                                                                     | PENTRA | PENTRA                           |                         |
|                          |                                                        | ES60<br>Care ST | CRP 200 |        |        | 60 8P - 60 18P<br>ADVIA 60 |                |                                                                                            |        |                                  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 7.9             | 7.8     |        |        | 8.0                        |                |                                                                                            |        |                                  | ± 0.8                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 4.60            | 4.55    |        |        | 4.62                       |                |                                                                                            |        |                                  | ± 0.18                  |
| HB HGB                   | g/dl                                                   | 13.5            | 13.5    |        |        | 13.6                       |                |                                                                                            |        |                                  | ± 0.5                   |
|                          | g/l                                                    | 135             | 135     |        |        | 136                        |                |                                                                                            |        |                                  | ± 5                     |
|                          | mmol/l                                                 | 8.38            | 8.38    |        |        | 8.45                       |                |                                                                                            |        |                                  | ± 0.31                  |
| HT HCT                   | %                                                      | 37.7            | 37.3    |        |        | 37.4                       |                |                                                                                            |        |                                  | ± 2.5                   |
|                          | l/l                                                    | 0.377           | 0.373   |        |        | 0.374                      |                |                                                                                            |        |                                  | ± 0.025                 |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 82              | 82      |        |        | 81                         |                |                                                                                            |        |                                  | ± 4                     |
| TGMH MCH                 | pg                                                     | 29.3            | 29.7    |        |        | 29.4                       |                |                                                                                            |        |                                  | ± 2.0                   |
|                          | fmol                                                   | 1.82            | 1.84    |        |        | 1.83                       |                |                                                                                            |        |                                  | ± 0.12                  |
| CCMH MCHC                | g/dl                                                   | 35.8            | 36.2    |        |        | 36.3                       |                |                                                                                            |        |                                  | ± 3.0                   |
|                          | g/l                                                    | 358             | 362     |        |        | 363                        |                |                                                                                            |        |                                  | ± 30                    |
|                          | mmol/l                                                 | 22.23           | 22.47   |        |        | 22.57                      |                |                                                                                            |        |                                  | ± 1.86                  |
| IDR RDW                  | %                                                      | 14.5            | 14.0    |        |        | 13.5                       |                |                                                                                            |        |                                  | ± 3.0                   |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 258             | 253     |        |        | 260                        |                |                                                                                            |        |                                  | ± 40                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 8.4             | 8.3     |        |        | 7.8                        |                |                                                                                            |        |                                  | ± 2.0                   |
| LYMPHO                   | #                                                      | 2.92            | 2.96    |        |        | 2.92                       |                |                                                                                            |        |                                  | ± 0.40                  |
|                          | %                                                      | 37.0            | 38.0    |        |        | 36.5                       |                |                                                                                            |        |                                  | ± 5.0                   |
| MONO                     | #                                                      | 0.47            | 0.55    |        |        | 0.48                       |                |                                                                                            |        |                                  | ± 0.30                  |
|                          | %                                                      | 6.0             | 7.0     |        |        | 6.0                        |                |                                                                                            |        |                                  | ± 4.0                   |
| GRANULO                  | #                                                      | 4.50            | 4.29    |        |        | 4.60                       |                |                                                                                            |        |                                  | ± 0.60                  |
|                          | %                                                      | 57.0            | 55.0    |        |        | 57.5                       |                |                                                                                            |        |                                  | ± 6.0                   |

Ref: TEMP-0830 Rev.41 FRONT / RECTO 130090307



| LOT                      |      | MX 430 H                                               |       | CONTROL  |  | H      |  |  (Exp.) |  | 2021-10-05<br>( YYYY - MM - DD ) |  |                         |  |                |  |        |  |        |         |
|--------------------------|------|--------------------------------------------------------|-------|----------|--|--------|--|--------------------------------------------------------------------------------------------|--|----------------------------------|--|-------------------------|--|----------------|--|--------|--|--------|---------|
|                          |      | Rev 1                                                  |       | ABX Lyse |  |        |  |                                                                                            |  |                                  |  |                         |  |                |  |        |  |        |         |
| PARAMETRES<br>PARAMETERS |      | UNITES<br>UNITS                                        |       |          |  |        |  |                                                                                            |  |                                  |  | TOLERANCES<br>TOLERANCE |  |                |  |        |  |        |         |
|                          |      |                                                        |       | MICROS   |  | MICROS |  | MICROS                                                                                     |  | MICROS                           |  |                         |  | ARGOS - HELIOS |  | PENTRA |  | PENTRA |         |
|                          |      | ES60<br>Care ST                                        |       | CRP 200  |  |        |  |                                                                                            |  | 60 8P - 60 18P<br>ADVIA 60       |  |                         |  |                |  |        |  |        |         |
| GB                       | WBC  | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 20.7  | 20.3     |  |        |  |                                                                                            |  | 21.0                             |  |                         |  |                |  |        |  |        | ± 1.6   |
| GR                       | RBC  | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 5.61  | 5.55     |  |        |  |                                                                                            |  | 5.64                             |  |                         |  |                |  |        |  |        | ± 0.20  |
| HB                       | HGB  | g/dl                                                   | 17.8  | 17.7     |  |        |  |                                                                                            |  | 18.0                             |  |                         |  |                |  |        |  |        | ± 0.6   |
|                          |      | g/l                                                    | 178   | 177      |  |        |  |                                                                                            |  | 180                              |  |                         |  |                |  |        |  |        | ± 6     |
|                          |      | mmol/l                                                 | 11.05 | 10.99    |  |        |  |                                                                                            |  | 11.18                            |  |                         |  |                |  |        |  |        | ± 0.37  |
| HT                       | HCT  | %                                                      | 49.9  | 49.4     |  |        |  |                                                                                            |  | 50.2                             |  |                         |  |                |  |        |  |        | ± 3.0   |
|                          |      | l/l                                                    | 0.499 | 0.494    |  |        |  |                                                                                            |  | 0.502                            |  |                         |  |                |  |        |  |        | ± 0.030 |
| VGM                      | MCV  | µm <sup>3</sup> ; fl                                   | 89    | 89       |  |        |  |                                                                                            |  | 89                               |  |                         |  |                |  |        |  |        | ± 4     |
| TGMH                     | MCH  | pg                                                     | 31.7  | 31.9     |  |        |  |                                                                                            |  | 31.9                             |  |                         |  |                |  |        |  |        | ± 2.0   |
|                          |      | fmol                                                   | 1.97  | 1.98     |  |        |  |                                                                                            |  | 1.98                             |  |                         |  |                |  |        |  |        | ± 0.12  |
| CCMH                     | MCHC | g/dl                                                   | 35.7  | 35.8     |  |        |  |                                                                                            |  | 35.9                             |  |                         |  |                |  |        |  |        | ± 3.0   |
|                          |      | g/l                                                    | 357   | 358      |  |        |  |                                                                                            |  | 359                              |  |                         |  |                |  |        |  |        | ± 30    |
|                          |      | mmol/l                                                 | 22.14 | 22.25    |  |        |  |                                                                                            |  | 22.27                            |  |                         |  |                |  |        |  |        | ± 1.86  |
| IDR                      | RDW  | %                                                      | 13.5  | 13.5     |  |        |  |                                                                                            |  | 13.0                             |  |                         |  |                |  |        |  |        | ± 3.0   |
| PLAQ.                    | PLTS | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 501   | 490      |  |        |  |                                                                                            |  | 500                              |  |                         |  |                |  |        |  |        | ± 45    |
| VPM                      | MPV  | µm <sup>3</sup> ; fl                                   | 8.0   | 8.0      |  |        |  |                                                                                            |  | 7.5                              |  |                         |  |                |  |        |  |        | ± 2.0   |
| LYMPHO                   | #    | 3.52                                                   | 3.45  |          |  |        |  |                                                                                            |  | 3.36                             |  |                         |  |                |  |        |  |        | ± 1.00  |
|                          | %    | 17.0                                                   | 17.0  |          |  |        |  |                                                                                            |  | 16.0                             |  |                         |  |                |  |        |  |        | ± 5.0   |
| MONO                     | #    | 1.04                                                   | 1.32  |          |  |        |  |                                                                                            |  | 0.95                             |  |                         |  |                |  |        |  |        | ± 0.80  |
|                          | %    | 5.0                                                    | 6.5   |          |  |        |  |                                                                                            |  | 4.5                              |  |                         |  |                |  |        |  |        | ± 4.0   |
| GRANULO                  | #    | 16.15                                                  | 15.53 |          |  |        |  |                                                                                            |  | 16.70                            |  |                         |  |                |  |        |  |        | ± 1.40  |
|                          | %    | 78.0                                                   | 76.5  |          |  |        |  |                                                                                            |  | 79.5                             |  |                         |  |                |  |        |  |        | ± 7.0   |