



LOT		MX 014 L		CONTROL		L		 (Exp.)		2014-04-05 (YYYY - MM - DD)											
		Rev 1																			
PARAMETRES PARAMETERS		UNITES UNITS		Minilysebio								TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l		2.0	2.1			2.0										± 0.4			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.32	2.35			2.34										± 0.15			
HB	HGB	g/dl		5.9	6.0			5.6										± 0.4			
		g/l		59	60			56										± 4			
		mmol/l		3.65	3.73			3.47										± 0.25			
HT	HCT	%		16.0	16.2			16.4										± 2.0			
		l/l		0.160	0.162			0.164										± 0.020			
VGM	MCV	µm ³ ; fl		69	69			70										± 4			
TGMH	MCH	pg		25.3	25.5			23.8										± 2.0			
		fmol		1.57	1.59			1.48										± 0.12			
CCMH	MCHC	g/dl		36.7	37.0			34.1										± 3.0			
		g/l		367	370			341										± 30			
		mmol/l		22.81	22.98			21.15										± 1.86			
IDR	RDW	%		13.7	13.7			13.6										± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		74	70			75										± 20			
VPM	MPV	µm ³ ; fl		8.3	8.2			7.7										± 2.0			
LYMPHO	#			1.31	1.3			1.27										± 0.40			
	%			65.7	64.0			63.5										± 7.0			
MONO	#			0.20	0.2			0.20										± 0.20			
	%			10.0	8.5			10.0										± 6.0			
GRANULO	#			0.49	0.6			0.53										± 0.40			
	%			24.5	27.5			26.5										± 7.0			



LOT		MX 014 N		Rev 1		CONTROL		N		 (Exp.)		2014-04-05 (YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
			CRP CRP 200	EMI CRP		60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		7.6	7.7		7.3						± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.60	4.59		4.64						± 0.18	
HB HGB	g/dl		13.4	13.4		13.2						± 0.5	
	g/l		134	134		132						± 5	
	mmol/l		8.32	8.32		8.22						± 0.31	
HT HCT	%		37.3	36.7		37.6						± 2.5	
	l/l		0.373	0.367		0.376						± 0.025	
VGM MCV	µm ³ ; fl		81	80		81						± 4	
TGMH MCH	pg		29.1	29.2		28.5						± 2.0	
	fmol		1.81	1.81		1.77						± 0.12	
CCMH MCHC	g/dl		36.0	36.5		35.2						± 3.0	
	g/l		360	365		352						± 30	
	mmol/l		22.33	22.66		21.86						± 1.86	
IDR RDW	%		13.0	13.0		13.8						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		237	235		244						± 40	
VPM MPV	µm ³ ; fl		7.9	8.3		7.5						± 2.0	
LYMPHO	#		2.63	2.6		2.41						± 0.40	
	%		34.6	33.5		33.2						± 5.0	
MONO	#		0.47	0.5		0.44						± 0.30	
	%		6.2	6.5		6.0						± 4.0	
GRANULO	#		4.50	4.6		4.44						± 0.60	
	%		59.2	60.0		61.0						± 6.0	



LOT

MX 014 H

Rev 1

CONTROL

H



2014-04-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	
			CRP CRP 200	EMI CRP		60 8P - 60 18P ADVIA 60				
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.3	20.5		19.5				± 2.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.66	5.67		5.73				± 0.20
HB HGB	g/dl		17.7	17.7		17.7				± 0.6
	g/l		177	177		177				± 6
	mmol/l		10.99	10.99		11.00				± 0.37
HT HCT	%		49.2	49.3		50.4				± 3.0
	l/l		0.492	0.493		0.504				± 0.030
VGM MCV	µm ³ ; fl		87	87		88				± 4
TGMH MCH	pg		31.3	31.2		30.9				± 2.0
	fmol		1.94	1.94		1.92				± 0.12
CCMH MCHC	g/dl		35.9	35.9		35.1				± 3.0
	g/l		359	359		351				± 30
	mmol/l		22.32	22.28		21.82				± 1.86
IDR RDW	%		12.7	12.3		13.2				± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		520	515		530				± 45
VPM MPV	µm ³ ; fl		7.7	8.3		7.4				± 2.0
LYMPHO	#		3.14	3.4		2.91				± 1.00
	%		15.5	16.5		14.9				± 5.0
MONO	#		1.32	1.1		1.07				± 0.80
	%		6.5	5.5		5.5				± 4.0
GRANULO	#		15.80	16.0		15.44				± 1.40
	%		77.8	78.0		79.2				± 7.0