



LOT		MX 034 L		CONTROL		L		 (Exp.)		2014-06-05 (YYYY - MM - DD)											
		Rev 1		Minilysebio																	
PARAMETRES PARAMETERS		UNITES UNITS										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
						CRP CRP 200		EMI CRP				60 8P - 60 18P ADVIA 60									
GB WBC		10 ³ /mm ³ ; 10 ⁹ /l				2.0		2.0				2.0								± 0.4	
GR RBC		10 ⁶ /mm ³ ; 10 ¹² /l				2.37		2.42				2.39								± 0.15	
HB HGB		g/dl				5.9		6.1				5.6								± 0.4	
		g/l				59		61				56								± 4	
		mmol/l				3.65		3.79				3.47								± 0.25	
HT HCT		%				16.1		16.5				16.5								± 2.0	
		l/l				0.161		0.165				0.165								± 0.020	
VGM MCV		µm ³ ; fl				68		68				69								± 4	
TGMH MCH		pg				24.8		25.2				23.3								± 2.0	
		fmol				1.54		1.57				1.45								± 0.12	
CCMH MCHC		g/dl				36.5		37.1				33.8								± 3.0	
		g/l				365		371				338								± 30	
		mmol/l				22.66		23.02				21.01								± 1.86	
IDR RDW		%				13.5		13.0				13.5								± 3.0	
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l				73		70				72								± 20	
VPM MPV		µm ³ ; fl				8.8		8.6				8.3								± 2.0	
LYMPHO		#				1.29		1.2				1.26								± 0.40	
		%				64.6		62.0				63.0								± 7.0	
MONO		#				0.20		0.2				0.20								± 0.20	
		%				10.0		9.0				10.0								± 6.0	
GRANULO		#				0.51		0.6				0.54								± 0.40	
		%				25.5		29.0				27.0								± 7.0	



LOT		MX 034 N		CONTROL		N		 (Exp.)		2014-06-05 (YYYY-MM-DD)			
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
			CRP CRP 200	EMI CRP		60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		7.1	7.5		7.1					± 0.8		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.53	4.54		4.58					± 0.18		
HB HGB	g/dl		13.3	13.3		13.1					± 0.5		
	g/l		133	133		131					± 5		
	mmol/l		8.26	8.26		8.15					± 0.31		
HT HCT	%		36.2	36.3		37.1					± 2.5		
	l/l		0.362	0.363		0.371					± 0.025		
VGM MCV	µm ³ ; fl		80	80		81					± 4		
TGMH MCH	pg		29.4	29.3		28.7					± 2.0		
	fmol		1.82	1.82		1.78					± 0.12		
CCMH MCHC	g/dl		36.7	36.6		35.4					± 3.0		
	g/l		367	366		354					± 30		
	mmol/l		22.79	22.74		21.98					± 1.86		
IDR RDW	%		13.0	12.5		13.0					± 3.0		
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		238	238		242					± 40		
VPM MPV	µm ³ ; fl		8.0	8.2		7.6					± 2.0		
LYMPHO	#		2.52	2.6		2.38					± 0.40		
	%		35.5	34.5		33.7					± 5.0		
MONO	#		0.47	0.5		0.46					± 0.30		
	%		6.7	6.5		6.5					± 4.0		
GRANULO	#		4.11	4.4		4.25					± 0.60		
	%		57.8	59.0		60.0					± 6.0		

Explore the future