

ABX Minotrol CRP



LOT

MC 402

CONTROL



(Exp.) 2017-04-05

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse																
				CONTROL				1	CONTROL				2	CONTROL				3	TOLERANCES TOLERANCE	
				MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE		
				CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0			± 0.4	7.6	7.6	7.6			± 1.0	20.1	20.1	20.1			± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.42	2.42	2.42			± 0.15	4.72	4.72	4.72			± 0.18	5.65	5.65	5.62			± 0.25
HB	HGB	g/dl	6.3	6.3	6.1			± 0.4	13.7	13.7	13.6			± 0.5	17.8	17.8	17.6			± 0.7
		g/l	63	63	61			± 4	137	137	136			± 5	178	178	176			± 7
		mmol/l	3.91	3.91	3.79			± 0.25	8.51	8.51	8.45			± 0.31	11.05	11.05	10.93			± 0.44
HT	HCT	%	16.9	16.9	16.9			± 2.0	37.8	37.8	37.8			± 2.5	50.3	50.3	49.5			± 3.0
		l/l	0.169	0.169	0.169			± 0.020	0.378	0.378	0.378			± 0.025	0.503	0.503	0.495			± 0.030
VGM	MCV	µm ³ ; fl	70	70	70			± 4	80	80	80			± 4	89	89	88			± 4
TGMH	MCH	pg	26.0	26.0	25.2			± 2.4	29.0	29.0	28.8			± 2.4	31.5	31.5	31.3			± 2.8
		fmol	1.62	1.62	1.57			± 0.15	1.80	1.80	1.79			± 0.15	1.96	1.96	1.94			± 0.18
CCMH	MCHC	g/dl	37.2	37.2	36.0			± 3.0	36.3	36.3	36.0			± 3.0	35.4	35.4	35.6			± 3.0
		g/l	372	372	360			± 30	363	363	360			± 30	354	354	356			± 30
		mmol/l	23.10	23.10	22.36			± 1.86	22.53	22.53	22.37			± 1.86	21.98	21.98	22.10			± 1.86
IDR	RDW	%	13.8	13.8	13.4			± 3.0	13.2	13.2	13.0			± 3.0	12.8	12.8	12.5			± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	70	70	75			± 20	255	255	250			± 40	510	510	500			± 60
VPM	MPV	µm ³ ; fl	9.0	9.0	8.0			± 2.0	8.7	8.7	8.0			± 2.0	8.6	8.6	7.8			± 2.0
LYMPHO		#	1.3	1.3	1.3			± 0.3	2.6	2.6	2.7			± 0.6	3.3	3.3	3.4			± 1.3
		%	64.0	64.0	64.0			± 8.0	34.0	34.0	35.0			± 6.0	16.5	16.5	17.0			± 6.0
MONO		#	0.2	0.2	0.2			± 0.2	0.5	0.5	0.5			± 0.4	1.1	1.1	1.0			± 1.0
		%	10.0	10.0	10.0			± 7.0	7.0	7.0	7.0			± 5.0	5.5	5.5	5.0			± 5.0
GRANULO		#	0.5	0.5	0.5			± 0.3	4.5	4.5	4.4			± 0.6	15.7	15.7	15.7			± 1.5
		%	26.0	26.0	26.0			± 8.0	59.0	59.0	58.0			± 6.0	78.0	78.0	78.0			± 7.0
CRP		mg/dl	0.69	0.54	N/A			± 0.30	2.48	2.27	N/A			± 0.70	4.49	4.28	N/A			± 1.00
		mg/l	6.9	5.4	N/A			± 3.0	24.8	22.7	N/A			± 7.0	44.9	42.8	N/A			± 10.0

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A