

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 403

CONTROL



(Exp.) 2017-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse																
				CONTROL				1	CONTROL				2	CONTROL				3	TOLERANCES TOLERANCE	
				MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE		
				CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0			± 0.4	7.7	7.7	7.7			± 1.0	20.3	20.3	20.6			± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.25	2.25	2.27			± 0.15	4.68	4.68	4.70			± 0.18	5.70	5.70	5.78			± 0.25
HB	HGB	g/dl	5.9	5.9	5.9			± 0.4	13.5	13.5	13.5			± 0.5	17.8	17.8	17.8			± 0.7
		g/l	59	59	59			± 4	135	135	135			± 5	178	178	178			± 7
		mmol/l	3.66	3.66	3.66			± 0.25	8.38	8.38	8.38			± 0.31	11.05	11.05	11.05			± 0.44
HT	HCT	%	16.0	16.0	15.9			± 2.0	37.9	37.9	37.6			± 2.5	50.7	50.7	50.9			± 3.0
		l/l	0.160	0.160	0.159			± 0.020	0.379	0.379	0.376			± 0.025	0.507	0.507	0.509			± 0.030
VGM	MCV	µm ³ ; fl	71	71	70			± 4	81	81	80			± 4	89	89	88			± 4
TGMH	MCH	pg	26.2	26.2	26.0			± 2.4	28.8	28.8	28.7			± 2.4	31.2	31.2	30.8			± 2.8
		fmol	1.63	1.63	1.61			± 0.15	1.79	1.79	1.78			± 0.15	1.94	1.94	1.91			± 0.18
CCMH	MCHC	g/dl	36.9	36.9	37.1			± 3.0	35.6	35.6	35.9			± 3.0	35.1	35.1	35.0			± 3.0
		g/l	369	369	371			± 30	356	356	359			± 30	351	351	350			± 30
		mmol/l	22.94	22.94	23.06			± 1.86	22.12	22.12	22.30			± 1.86	21.79	21.79	21.73			± 1.86
IDR	RDW	%	14.0	14.0	12.5			± 3.0	14.0	14.0	12.5			± 3.0	13.5	13.5	12.0			± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	265	265	266			± 40	515	515	520			± 60
VPM	MPV	µm ³ ; fl	9.2	9.2	8.6			± 2.0	8.7	8.7	8.2			± 2.0	8.4	8.4	7.9			± 2.0
LYMPHO		#	1.3	1.3	1.2			± 0.3	2.7	2.7	2.8			± 0.6	3.2	3.2	3.4			± 1.3
		%	62.5	62.5	62.0			± 7.0	34.5	34.5	36.0			± 6.0	16.0	16.0	16.5			± 6.0
MONO		#	0.2	0.2	0.2			± 0.2	0.7	0.7	0.5			± 0.4	1.2	1.2	1.0			± 1.0
		%	9.5	9.5	9.0			± 6.0	8.5	8.5	7.0			± 5.0	6.0	6.0	5.0			± 5.0
GRANULO		#	0.6	0.6	0.6			± 0.3	4.4	4.4	4.4			± 0.6	15.8	15.8	16.2			± 1.5
		%	28.0	28.0	29.0			± 8.0	57.0	57.0	57.0			± 6.0	78.0	78.0	78.5			± 7.0
CRP		mg/dl	0.95	0.76	N/A			± 0.30	2.73	2.49	N/A			± 0.70	4.35	3.94	N/A			± 1.00
		mg/l	9.5	7.6	N/A			± 3.0	27.3	24.9	N/A			± 7.0	43.5	39.4	N/A			± 10.0

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A