

# ABX Minotrol CRP



**LOT**

**MC 404**

Rev 2

**CONTROL**



(Exp.) 2017-10-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse																
		CONTROL				1	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE		MICROS				TOLERANCES TOLERANCE	
		CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P				
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.0	2.0	2.0		± 0.4	7.4	7.4	7.5		± 1.0	19.8	19.8	20.0		± 2.2		
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.43	2.43	2.44		± 0.15	4.77	4.77	4.76		± 0.18	5.76	5.76	5.79		± 0.25		
HB HGB	g/dl	6.2	6.2	6.3		± 0.4	13.4	13.4	13.5		± 0.5	17.9	17.9	18.0		± 0.7		
	g/l	62	62	63		± 4	134	134	135		± 5	179	179	180		± 7		
HT HCT	mmol/l	3.85	3.85	3.91		± 0.25	8.32	8.32	8.38		± 0.31	11.12	11.12	11.18		± 0.44		
	%	17.0	17.0	17.1		± 2.0	38.2	38.2	38.1		± 2.5	51.3	51.3	51.5		± 3.0		
VGM MCV	l/l	0.170	0.170	0.171		± 0.020	0.382	0.382	0.381		± 0.025	0.513	0.513	0.515		± 0.030		
	µm <sup>3</sup> ; fl	70	70	70		± 4	80	80	80		± 4	89	89	89		± 4		
TGMH MCH	pg	25.5	25.5	25.8		± 2.4	28.1	28.1	28.4		± 2.4	31.1	31.1	31.1		± 2.8		
	fmol	1.58	1.58	1.60		± 0.15	1.74	1.74	1.76		± 0.15	1.93	1.93	1.93		± 0.18		
CCMH MCHC	g/dl	36.4	36.4	36.9		± 3.0	35.1	35.1	35.5		± 3.0	34.9	34.9	34.9		± 3.0		
	g/l	364	364	369		± 30	351	351	355		± 30	349	349	349		± 30		
IDR RDW	mmol/l	22.63	22.63	22.91		± 1.86	21.81	21.81	22.02		± 1.86	21.68	21.68	21.69		± 1.86		
	%	14.0	14.0	12.6		± 3.0	14.5	14.5	13.4		± 3.0	13.5	13.5	12.0		± 3.0		
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	75	75	76		± 20	250	250	252		± 40	485	485	490		± 60		
VPM MPV	µm <sup>3</sup> ; fl	9.0	9.0	7.8		± 2.0	8.3	8.3	7.5		± 2.0	7.6	7.6	6.8		± 2.0		
LYMPHO	#	1.3	1.3	1.3		± 0.2	2.5	2.5	2.6		± 0.5	3.1	3.1	3.2		± 1.2		
	%	65.5	65.5	64.0		± 7.0	33.5	33.5	34.0		± 6.0	15.5	15.5	16.0		± 6.0		
MONO	#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.5		± 0.4	1.1	1.1	0.9		± 0.9		
	%	8.5	8.5	8.5		± 6.0	8.5	8.5	7.0		± 5.0	5.5	5.5	4.5		± 4.5		
GRANULO	#	0.5	0.5	0.6		± 0.2	4.3	4.3	4.4		± 0.5	15.6	15.6	15.9		± 1.4		
	%	26.0	26.0	27.5		± 8.0	58.0	58.0	59.0		± 6.0	79.0	79.0	79.5		± 7.0		
CRP	mg/dl	0.76	0.54	N/A		± 0.30	2.70	2.32	N/A		± 0.70	4.99	4.78	N/A		± 1.00		
	mg/l	7.6	5.4	N/A		± 3.0	27.0	23.2	N/A		± 7.0	49.9	47.8	N/A		± 10.0		

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A