

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MC 405**

**CONTROL**



**(Exp.) 2018-01-05**  
(YYYY - MM - DD)

Rev 1

		<b>ABX Lyse</b>																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS						MICROS						MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.1	2.1	2.2		± 0.4	7.6	7.6	7.6		± 1.0	20.2	20.2	20.3		± 2.2			
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.21	2.21	2.23		± 0.15	4.67	4.67	4.67		± 0.18	5.64	5.64	5.69		± 0.25			
HB HGB	g/dl	5.8	5.8	5.8		± 0.4	13.5	13.5	13.6		± 0.5	17.7	17.7	17.8		± 0.7			
	g/l	58	58	58		± 4	135	135	136		± 5	177	177	178		± 7			
HT HCT	mmol/l	3.60	3.60	3.60		± 0.25	8.38	8.38	8.45		± 0.31	10.99	10.99	11.05		± 0.44			
	%	15.5	15.5	15.6		± 2.0	37.8	37.8	37.8		± 2.5	49.6	49.6	50.6		± 3.0			
VGM MCV	l/l	0.155	0.155	0.156		± 0.020	0.378	0.378	0.378		± 0.025	0.496	0.496	0.506		± 0.030			
	µm <sup>3</sup> ; fl	70	70	70		± 4	81	81	81		± 4	88	88	89		± 4			
TGMH MCH	pg	26.2	26.2	26.0		± 2.4	28.9	28.9	29.1		± 2.4	31.4	31.4	31.3		± 2.8			
	fmol	1.63	1.63	1.62		± 0.15	1.80	1.80	1.81		± 0.15	1.95	1.95	1.94		± 0.18			
CCMH MCHC	g/dl	37.5	37.5	37.2		± 3.0	35.7	35.7	36.0		± 3.0	35.7	35.7	35.1		± 3.0			
	g/l	375	375	372		± 30	357	357	360		± 30	357	357	351		± 30			
	mmol/l	23.28	23.28	23.07		± 1.86	22.16	22.16	22.33		± 1.86	22.15	22.15	21.83		± 1.86			
IDR RDW	%	14.5	14.5	13.0		± 3.0	14.0	14.0	13.5		± 3.0	13.5	13.5	12.5		± 3.0			
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	69	69	71		± 20	255	255	257		± 40	530	530	535		± 60			
VPM MPV	µm <sup>3</sup> ; fl	8.7	8.7	8.0		± 2.0	8.1	8.1	7.7		± 2.0	8.3	8.3	7.8		± 2.0			
LYMPHO	#	1.4	1.4	1.5		± 0.2	2.7	2.7	2.7		± 0.5	3.2	3.2	3.2		± 1.3			
	%	67.5	67.5	66.0		± 8.0	35.0	35.0	35.0		± 6.0	16.0	16.0	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.5	0.5	0.5		± 0.4	1.1	1.1	1.0		± 1.0			
	%	8.5	8.5	8.5		± 7.0	7.0	7.0	6.5		± 5.0	5.5	5.5	5.0		± 5.0			
GRANULO	#	0.6	0.6	0.6		± 0.2	4.4	4.4	4.4		± 0.5	15.9	15.9	16.0		± 1.5			
	%	24.0	24.0	25.5		± 7.0	58.0	58.0	58.5		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.61	0.46	N/A		± 0.30	2.24	1.94	N/A		± 0.70	3.69	3.27	N/A		± 1.00			
	mg/l	6.1	4.6	N/A		± 3.0	22.4	19.4	N/A		± 7.0	36.9	32.7	N/A		± 10.0			

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A