

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 043

CONTROL



(Exp.) 2013-07-05

Rev 2

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse													TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3			
				MICROS				MICROS				MICROS						
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0		± 0.4	7.6	7.6	7.7		± 1.0	20.5	20.5	20.6		± 2.2	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.35	2.33		± 0.15	4.60	4.60	4.58		± 0.18	5.72	5.72	5.72		± 0.25	
HB	HGB	g/dl	6.1	6.1	6.1		± 0.4	13.5	13.5	13.5		± 0.5	18.0	18.0	18.0		± 0.7	
		g/l	61	61	61		± 4	135	135	135		± 5	180	180	180		± 7	
		mmol/l	3.79	3.79	3.79		± 0.25	8.38	8.38	8.38		± 0.31	11.18	11.18	11.18		± 0.44	
HT	HCT	%	16.7	16.7	16.5		± 2.0	37.7	37.7	37.6		± 2.5	50.9	50.9	50.9		± 3.0	
		l/l	0.167	0.167	0.165		± 0.020	0.377	0.377	0.376		± 0.025	0.509	0.509	0.509		± 0.030	
VGM	MCV	µm ³ ; fl	71	71	71		± 4	82	82	82		± 4	89	89	89		± 4	
TGMH	MCH	pg	26.0	26.0	26.2		± 2.4	29.3	29.3	29.5		± 2.4	31.5	31.5	31.5		± 2.8	
		fmol	1.61	1.61	1.63		± 0.15	1.82	1.82	1.83		± 0.15	1.95	1.95	1.95		± 0.18	
CCMH	MCHC	g/dl	36.6	36.6	36.9		± 3.0	35.8	35.8	35.9		± 3.0	35.4	35.4	35.4		± 3.0	
		g/l	366	366	369		± 30	358	358	359		± 30	354	354	354		± 30	
		mmol/l	22.70	22.70	22.90		± 1.86	22.23	22.23	22.32		± 1.86	21.96	21.96	21.96		± 1.86	
IDR	RDW	%	13.6	13.6	13.8		± 3.0	13.8	13.8	13.2		± 3.0	13.0	13.0	13.0		± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	73	73	74		± 20	260	260	260		± 40	515	515	510		± 60	
VPM	MPV	µm ³ ; fl	8.5	8.5	8.3		± 2.0	7.8	7.8	7.5		± 2.0	7.8	7.8	7.6		± 2.0	
LYMPHO		#	1.3	1.3	1.3		± 0.3	2.5	2.5	2.5		± 0.5	3.0	3.0	3.0		± 1.1	
		%	64.0	64.0	62.5		± 7.0	32.3	32.3	32.0		± 6.0	14.5	14.5	14.5		± 5.0	
MONO		#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6		± 0.4	1.1	1.1	1.0		± 1.0	
		%	9.5	9.5	9.5		± 7.0	8.4	8.4	8.0		± 5.0	5.5	5.5	5.0		± 5.0	
GRANULO		#	0.5	0.5	0.6		± 0.3	4.5	4.5	4.6		± 0.5	16.4	16.4	16.6		± 1.5	
		%	26.5	26.5	28.0		± 7.0	59.3	59.3	60.0		± 6.0	80.0	80.0	80.5		± 7.0	
CRP		mg/dl	1.00	0.60	N/A		± 0.30	2.80	2.10	N/A		± 0.70	4.80	4.00	N/A		± 1.00	
		mg/l	10.0	6.0	N/A		± 3.0	28.0	21.0	N/A		± 7.0	48.0	40.0	N/A		± 10.0	

9930074-A
FRONT/RECTO
Ref: TEMP-1070 Rev.31