

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 016

CONTROL

(Exp.) **2016-04-05**
(YYYY - MM - DD)

Rev 3

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					1	CONTROL					2	CONTROL					3
		MICROS					TOLERANCES TOLERANCE	MICROS					TOLERANCES TOLERANCE	MICROS					TOLERANCES TOLERANCE
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0			± 0.4	7.2	7.2	7.3			± 1.0	19.3	19.3	19.7			± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.45	2.45	2.44			± 0.15	4.65	4.65	4.65			± 0.18	5.62	5.62	5.62			± 0.25
HB HGB	g/dl	6.1	6.1	6.1			± 0.4	13.4	13.4	13.5			± 0.5	18.0	18.0	18.0			± 0.7
	g/l	61	61	61			± 4	134	134	135			± 5	180	180	180			± 7
HT HCT	mmol/l	3.79	3.79	3.79			± 0.25	8.32	8.32	8.38			± 0.31	11.18	11.18	11.18			± 0.44
	%	16.9	16.9	16.6			± 2.0	37.2	37.2	37.2			± 2.5	50.0	50.0	50.0			± 3.0
VGM MCV	l/l	0.169	0.169	0.166			± 0.020	0.372	0.372	0.372			± 0.025	0.500	0.500	0.500			± 0.030
	µm ³ ; fl	69	69	68			± 4	80	80	80			± 4	89	89	89			± 4
TGMH MCH	pg	24.9	24.9	25.0			± 2.4	28.8	28.8	29.0			± 2.4	32.0	32.0	32.0			± 2.8
	fmol	1.55	1.55	1.55			± 0.15	1.79	1.79	1.80			± 0.15	1.99	1.99	1.99			± 0.18
CCMH MCHC	g/dl	36.1	36.1	36.8			± 3.0	36.0	36.0	36.3			± 3.0	36.0	36.0	36.0			± 3.0
	g/l	361	361	368			± 30	360	360	363			± 30	360	360	360			± 30
	mmol/l	22.41	22.41	22.83			± 1.86	22.37	22.37	22.54			± 1.86	22.35	22.35	22.35			± 1.86
IDR RDW	%	13.5	13.5	14.0			± 3.0	13.0	13.0	13.5			± 3.0	12.5	12.5	13.0			± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	250	250	252			± 40	495	495	496			± 60
VPM MPV	µm ³ ; fl	9.0	9.0	8.5			± 2.0	8.0	8.0	7.5			± 2.0	7.8	7.8	7.5			± 2.0
LYMPHO	#	1.3	1.3	1.3			± 0.4	2.4	2.4	2.5			± 0.5	2.9	2.9	3.1			± 1.2
	%	67.0	67.0	65.0			± 8.0	34.0	34.0	33.5			± 6.0	15.0	15.0	15.5			± 6.0
MONO	#	0.2	0.2	0.2			± 0.2	0.6	0.6	0.6			± 0.4	1.2	1.2	1.1			± 1.0
	%	10.5	10.5	11.0			± 7.0	8.5	8.5	8.5			± 5.0	6.0	6.0	5.5			± 5.0
GRANULO	#	0.5	0.5	0.5			± 0.4	4.2	4.2	4.2			± 0.5	15.2	15.2	15.5			± 1.4
	%	22.5	22.5	24.0			± 8.0	57.5	57.5	58.0			± 6.0	79.0	79.0	79.0			± 7.0
CRP	mg/dl	0.71	0.60	N/A			± 0.30	2.63	2.53	N/A			± 0.70	4.32	4.20	N/A			± 1.00
	mg/l	7.1	6.0	N/A			± 3.0	26.3	25.3	N/A			± 7.0	43.2	42.0	N/A			± 10.0

Ref: TEMP-1070 Rev.36 FRONT/RECTO 9930089-A