

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 073

CONTROL



(Exp.) 2013-10-05

(YYYY - MM - DD)

Rev 1

		ABX Lyse																				
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					1	TOLERANCES TOLERANCE	CONTROL					2	TOLERANCES TOLERANCE	CONTROL					3	TOLERANCES TOLERANCE
		MICROS							MICROS							MICROS						
		CRP	CRP 200	60 8P 60 18P					CRP	CRP 200	60 8P 60 18P					CRP	CRP 200	60 8P 60 18P				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	1.9	1.9	2.0			± 0.4	7.3	7.3	7.6			± 1.0	19.5	19.5	19.7			± 2.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.35	2.36			± 0.15	4.58	4.58	4.58			± 0.18	5.51	5.51	5.54			± 0.25			
HB HGB	g/dl	6.0	6.0	6.0			± 0.4	13.2	13.2	13.2			± 0.5	17.7	17.7	17.8			± 0.7			
	g/l	60	60	60			± 4	132	132	132			± 5	177	177	178			± 7			
	mmol/l	3.73	3.73	3.73			± 0.25	8.20	8.20	8.20			± 0.31	10.99	10.99	11.05			± 0.44			
HT HCT	%	16.5	16.5	16.5			± 2.0	37.1	37.1	37.1			± 2.5	49.6	49.6	49.9			± 3.0			
	l/l	0.165	0.165	0.165			± 0.020	0.371	0.371	0.371			± 0.025	0.496	0.496	0.499			± 0.030			
VGM MCV	µm ³ , fl	70	70	70			± 4	81	81	81			± 4	90	90	90			± 4			
TGMH MCH	pg	25.5	25.5	25.4			± 2.4	28.8	28.8	28.8			± 2.4	32.1	32.1	32.1			± 2.8			
	fmol	1.59	1.59	1.58			± 0.15	1.79	1.79	1.79			± 0.15	1.99	1.99	2.00			± 0.18			
CCMH MCHC	g/dl	36.5	36.5	36.3			± 3.0	35.6	35.6	35.6			± 3.0	35.7	35.7	35.7			± 3.0			
	g/l	365	365	363			± 30	356	356	356			± 30	357	357	357			± 30			
	mmol/l	22.65	22.65	22.55			± 1.86	22.10	22.10	22.10			± 1.86	22.17	22.17	22.17			± 1.86			
IDR RDW	%	13.8	13.8	13.9			± 3.0	13.6	13.6	13.4			± 3.0	12.7	12.7	12.8			± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	250	250	250			± 40	500	500	500			± 60			
VPM MPV	µm ³ , fl	9.1	9.1	8.9			± 2.0	8.1	8.1	7.8			± 2.0	8.1	8.1	7.9			± 2.0			
LYMPHO	#	1.2	1.2	1.3			± 0.3	2.4	2.4	2.6			± 0.5	2.8	2.8	3.0			± 1.2			
	%	64.0	64.0	62.5			± 8.0	33.5	33.5	34.5			± 6.0	14.5	14.5	15.0			± 6.0			
MONO	#	0.2	0.2	0.2			± 0.2	0.5	0.5	0.5			± 0.4	1.1	1.1	1.0			± 1.0			
	%	10.0	10.0	10.5			± 7.0	7.5	7.5	6.5			± 5.0	5.5	5.5	5.0			± 5.0			
GRANULO	#	0.5	0.5	0.5			± 0.2	4.3	4.3	4.5			± 0.6	15.6	15.6	15.8			± 1.4			
	%	26.0	26.0	27.0			± 7.0	59.0	59.0	59.0			± 6.0	80.0	80.0	80.0			± 7.0			
CRP	mg/dl	1.08	0.69	N/A			± 0.30	2.85	2.24	N/A			± 0.70	4.77	4.03	N/A			± 1.00			
	mg/l	10.8	6.9	N/A			± 3.0	28.5	22.4	N/A			± 7.0	47.7	40.3	N/A			± 10.0			

9930074-A FRONT/RECTO Ref: TEMP-1070 Rev.31