

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 044

CONTROL



(Exp.) 2014-07-05

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse													TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3			
				MICROS				MICROS				MICROS						
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0		± 0.4	7.4	7.4	7.5		± 1.0	20.1	20.1	20.5		± 2.2	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.37	2.37	2.39		± 0.15	4.60	4.60	4.64		± 0.18	5.66	5.66	5.72		± 0.25	
HB	HGB	g/dl	6.1	6.1	6.0		± 0.4	13.5	13.5	13.5		± 0.5	17.9	17.9	17.9		± 0.7	
		g/l	61	61	60		± 4	135	135	135		± 5	179	179	179		± 7	
		mmol/l	3.79	3.79	3.73		± 0.25	8.38	8.38	8.38		± 0.31	11.12	11.12	11.12		± 0.44	
HT	HCT	%	16.4	16.4	16.5		± 2.0	37.3	37.3	37.6		± 2.5	50.4	50.4	50.9		± 3.0	
		l/l	0.164	0.164	0.165		± 0.020	0.373	0.373	0.376		± 0.025	0.504	0.504	0.509		± 0.030	
VGM	MCV	µm ³ ; fl	69	69	69		± 4	81	81	81		± 4	89	89	89		± 4	
TGMH	MCH	pg	25.7	25.7	25.1		± 2.4	29.3	29.3	29.1		± 2.4	31.6	31.6	31.3		± 2.8	
		fmol	1.60	1.60	1.56		± 0.15	1.82	1.82	1.81		± 0.15	1.96	1.96	1.94		± 0.18	
CCMH	MCHC	g/dl	37.3	37.3	36.4		± 3.0	36.2	36.2	35.9		± 3.0	35.5	35.5	35.2		± 3.0	
		g/l	373	373	364		± 30	362	362	359		± 30	355	355	352		± 30	
		mmol/l	23.16	23.16	22.59		± 1.86	22.50	22.50	22.31		± 1.86	22.07	22.07	21.84		± 1.86	
IDR	RDW	%	13.5	13.5	13.5		± 3.0	13.4	13.4	13.6		± 3.0	13.2	13.2	13.5		± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	79	79	79		± 20	255	255	255		± 40	500	500	505		± 60	
VPM	MPV	µm ³ ; fl	8.8	8.8	8.4		± 2.0	8.1	8.1	7.7		± 2.0	7.8	7.8	7.5		± 2.0	
LYMPHO		#	1.3	1.3	1.3		± 0.3	2.5	2.5	2.6		± 0.5	3.0	3.0	3.1		± 1.3	
		%	65.5	65.5	64.5		± 8.0	34.0	34.0	34.0		± 6.0	15.0	15.0	15.0		± 6.0	
MONO		#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6		± 0.4	1.4	1.4	1.1		± 1.0	
		%	9.5	9.5	9.0		± 7.0	8.5	8.5	7.5		± 5.0	7.0	7.0	5.5		± 5.0	
GRANULO		#	0.5	0.5	0.5		± 0.2	4.3	4.3	4.4		± 0.5	15.7	15.7	16.3		± 1.5	
		%	25.0	25.0	26.5		± 7.0	57.5	57.5	58.5		± 6.0	78.0	78.0	79.5		± 7.0	
CRP		mg/dl	0.90	0.50	N/A		± 0.30	2.90	2.50	N/A		± 0.70	4.60	4.10	N/A		± 1.00	
		mg/l	9.0	5.0	N/A		± 3.0	29.0	25.0	N/A		± 7.0	46.0	41.0	N/A		± 10.0	

9930074-A FRONT/RECTO Ref: TEMP-1070 Rev.34