

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 104

CONTROL



(Exp.) 2015-01-05

(YYYY - MM - DD)

Rev 1

		ABX Lyse																				
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					1	TOLERANCES TOLERANCE	CONTROL					2	TOLERANCES TOLERANCE	CONTROL					3	TOLERANCES TOLERANCE
		MICROS							MICROS							MICROS						
		CRP	CRP 200	60 8P 60 18P					CRP	CRP 200	60 8P 60 18P					CRP	CRP 200	60 8P 60 18P				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1	2.1			± 0.4	7.6	7.6	7.6			± 1.0	20.4	20.4	20.4			± 2.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.21	2.21	2.22			± 0.15	4.45	4.45	4.50			± 0.18	5.52	5.52	5.60			± 0.25			
HB HGB	g/dl	6.0	6.0	5.9			± 0.4	13.3	13.3	13.3			± 0.5	17.7	17.7	17.8			± 0.7			
	g/l	60	60	59			± 4	133	133	133			± 5	177	177	178			± 7			
	mmol/l	3.73	3.73	3.66			± 0.25	8.26	8.26	8.26			± 0.31	10.99	10.99	11.05			± 0.44			
HT HCT	%	15.9	15.9	15.8			± 2.0	36.5	36.5	36.9			± 2.5	49.7	49.7	50.4			± 3.0			
	l/l	0.159	0.159	0.158			± 0.020	0.365	0.365	0.369			± 0.025	0.497	0.497	0.504			± 0.030			
VGM MCV	µm ³ , fl	72	72	71			± 4	82	82	82			± 4	90	90	90			± 4			
TGMH MCH	pg	27.1	27.1	26.6			± 2.4	29.9	29.9	29.6			± 2.4	32.1	32.1	31.8			± 2.8			
	fmol	1.69	1.69	1.65			± 0.15	1.86	1.86	1.84			± 0.15	1.99	1.99	1.97			± 0.18			
CCMH MCHC	g/dl	37.7	37.7	37.4			± 3.0	36.4	36.4	36.0			± 3.0	35.6	35.6	35.3			± 3.0			
	g/l	377	377	374			± 30	364	364	360			± 30	356	356	353			± 30			
	mmol/l	23.42	23.42	23.25			± 1.86	22.63	22.63	22.38			± 1.86	22.13	22.13	21.93			± 1.86			
IDR RDW	%	13.7	13.7	14.0			± 3.0	13.5	13.5	13.5			± 3.0	13.0	13.0	13.0			± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	255	255	253			± 40	505	505	504			± 60			
VPM MPV	µm ³ , fl	8.3	8.3	7.8			± 2.0	7.7	7.7	7.5			± 2.0	7.7	7.7	7.5			± 2.0			
LYMPHO	#	1.4	1.4	1.3			± 0.2	2.6	2.6	2.6			± 0.5	2.9	2.9	2.9			± 1.3			
	%	64.5	64.5	64.0			± 8.0	34.0	34.0	34.0			± 6.0	14.0	14.0	14.0			± 6.0			
MONO	#	0.2	0.2	0.2			± 0.2	0.6	0.6	0.6			± 0.4	1.1	1.1	1.0			± 1.0			
	%	10.0	10.0	9.5			± 7.0	8.0	8.0	7.5			± 5.0	5.5	5.5	5.0			± 5.0			
GRANULO	#	0.5	0.5	0.6			± 0.2	4.4	4.4	4.4			± 0.5	16.4	16.4	16.5			± 1.5			
	%	25.5	25.5	26.5			± 7.0	58.0	58.0	58.5			± 6.0	80.5	80.5	81.0			± 7.0			
CRP	mg/dl	0.72	0.50	N/A			± 0.30	2.68	2.32	N/A			± 0.70	4.35	3.85	N/A			± 1.00			
	mg/l	7.2	5.0	N/A			± 3.0	26.8	23.2	N/A			± 7.0	43.5	38.5	N/A			± 10.0			

9930089-A FRONT/RECTO Ref: TEMP-1070 Rev.36