

ABX Minotrol CRP



LOT

MC 407

CONTROL



(Exp.) 2018-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		Minilysebio												TOLERANCES TOLERANCE							
				CONTROL				1		CONTROL				2				CONTROL				3	
				MICROS				TOLERANCES TOLERANCE		MICROS				TOLERANCES TOLERANCE				MICROS				TOLERANCES TOLERANCE	
				CRP 200	60 8P 60 18P					CRP 200	60 8P 60 18P							CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l		2.0	2.0			± 0.4		7.5	7.7			± 1.0		19.7	20.4			± 2.2			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l		2.32	2.31			± 0.15		4.62	4.62			± 0.18		5.83	5.86			± 0.25			
HB	HGB	g/dl		5.8	5.9			± 0.4		13.0	13.4			± 0.5		17.7	18.2			± 0.7			
		g/l		58	59			± 4		130	134			± 5		177	182			± 7			
HT	HCT	%		16.2	16.2			± 2.0		37.4	37.4			± 2.5		50.1	50.4			± 3.0			
		l/l		0.162	0.162			± 0.020		0.374	0.374			± 0.025		0.501	0.504			± 0.030			
VGM	MCV	µm ³ ; fl		70	70			± 4		81	81			± 4		86	86			± 4			
TGMH	MCH	pg		25.0	25.5			± 2.4		28.2	29.0			± 2.4		30.4	31.1			± 2.8			
		fmol		1.55	1.59			± 0.15		1.75	1.80			± 0.15		1.89	1.93			± 0.18			
CCMH	MCHC	g/dl		35.7	36.5			± 3.0		34.8	35.8			± 3.0		35.3	36.1			± 3.0			
		g/l		357	365			± 30		348	358			± 30		353	361			± 30			
		mmol/l		22.16	22.66			± 1.86		21.63	22.24			± 1.86		21.92	22.43			± 1.86			
IDR	RDW	%		14.0	13.8			± 3.0		14.0	14.2			± 3.0		13.0	13.4			± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l		79	86			± 20		256	265			± 40		509	520			± 60			
VPM	MPV	µm ³ ; fl		8.1	8.2			± 2.0		8.5	8.9			± 2.0		8.0	8.4			± 2.0			
LYMPHO		#		1.4	1.3			± 0.2		2.6	2.7			± 0.5		3.2	3.2			± 1.2			
		%		68.6	67.0			± 8.0		35.2	34.5			± 6.0		16.1	15.5			± 6.0			
MONO		#		0.2	0.2			± 0.2		0.5	0.4			± 0.4		0.9	0.8			± 0.8			
		%		8.8	9.0			± 7.0		6.7	5.5			± 5.0		4.5	4.0			± 4.0			
GRANULO		#		0.5	0.5			± 0.2		4.4	4.6			± 0.5		15.7	16.4			± 1.5			
		%		22.5	24.0			± 7.0		58.3	60.0			± 6.0		79.7	80.5			± 7.0			
CRP		mg/dl		0.50	N/A			± 0.30		2.65	N/A			± 0.70		3.88	N/A			± 1.00			
		mg/l		5.0	N/A			± 3.0		26.5	N/A			± 7.0		38.8	N/A			± 10.0			

1300027267-A
Ref: TEMP-1070 Rev.37