

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 400

CONTROL

(Exp.) **2016-10-05**
(YYYY - MM - DD)

Rev 1

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS						MICROS						MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0		± 0.4	7.3	7.3	7.5		± 1.0	19.5	19.5	19.8		± 2.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.35	2.36		± 0.15	4.60	4.60	4.63		± 0.18	5.76	5.76	5.80		± 0.25			
HB HGB	g/dl	6.2	6.2	6.1		± 0.4	13.6	13.6	13.6		± 0.5	18.1	18.1	18.1		± 0.7			
	g/l	62	62	61		± 4	136	136	136		± 5	181	181	181		± 7			
HT HCT	mmol/l	3.85	3.85	3.79		± 0.25	8.45	8.45	8.45		± 0.31	11.24	11.24	11.24		± 0.44			
	%	16.5	16.5	16.5		± 2.0	37.3	37.3	37.5		± 2.5	50.1	50.1	50.5		± 3.0			
VGM MCV	l/l	0.165	0.165	0.165		± 0.020	0.373	0.373	0.375		± 0.025	0.501	0.501	0.505		± 0.030			
	µm ³ ; fl	70	70	70		± 4	81	81	81		± 4	87	87	87		± 4			
TGMH MCH	pg	26.4	26.4	25.8		± 2.4	29.6	29.6	29.4		± 2.4	31.4	31.4	31.2		± 2.8			
	fmol	1.64	1.64	1.61		± 0.15	1.84	1.84	1.82		± 0.15	1.95	1.95	1.94		± 0.18			
CCMH MCHC	g/dl	37.7	37.7	36.9		± 3.0	36.5	36.5	36.3		± 3.0	36.1	36.1	35.9		± 3.0			
	g/l	377	377	369		± 30	365	365	363		± 30	361	361	359		± 30			
	mmol/l	23.41	23.41	22.93		± 1.86	22.67	22.67	22.52		± 1.86	22.43	22.43	22.28		± 1.86			
IDR RDW	%	13.6	13.6	13.4		± 3.0	13.2	13.2	13.3		± 3.0	12.8	12.8	13.0		± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	72		± 20	254	254	250		± 40	510	510	505		± 60			
VPM MPV	µm ³ ; fl	8.7	8.7	8.3		± 2.0	8.0	8.0	7.6		± 2.0	7.9	7.9	7.5		± 2.0			
LYMPHO	#	1.3	1.3	1.3		± 0.3	2.5	2.5	2.6		± 0.5	3.0	3.0	3.2		± 1.2			
	%	64.0	64.0	62.8		± 8.0	34.0	34.0	35.0		± 6.0	15.5	15.5	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.5		± 0.4	1.2	1.2	1.0		± 1.0			
	%	10.0	10.0	10.0		± 7.0	8.0	8.0	7.0		± 5.0	6.0	6.0	5.0		± 5.0			
GRANULO	#	0.5	0.5	0.5		± 0.3	4.2	4.2	4.4		± 0.5	15.3	15.3	15.6		± 1.4			
	%	26.0	26.0	27.2		± 8.0	58.0	58.0	58.0		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.54	0.49	N/A		± 0.30	2.65	2.65	N/A		± 0.70	4.73	4.75	N/A		± 1.00			
	mg/l	5.4	4.9	N/A		± 3.0	26.5	26.5	N/A		± 7.0	47.3	47.5	N/A		± 10.0			

Ref: TEMP-1070 Rev.36 FRONT/RECTO 9930089-A