

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 104

CONTROL



(Exp.) **2015-01-05**

(YYYY - MM - DD)

Rev 1

		ABX Lyse																		
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					1	CONTROL					2	CONTROL					3	TOLERANCES TOLERANCE
		MICROS					TOLERANCES TOLERANCE	MICROS					TOLERANCES TOLERANCE	MICROS					TOLERANCES TOLERANCE	
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1	2.1			± 0.4	7.6	7.6	7.6			± 1.0	20.4	20.4	20.4			± 2.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.21	2.21	2.22			± 0.15	4.45	4.45	4.50			± 0.18	5.52	5.52	5.60			± 0.25	
HB HGB	g/dl	6.0	6.0	5.9			± 0.4	13.3	13.3	13.3			± 0.5	17.7	17.7	17.8			± 0.7	
	g/l	60	60	59			± 4	133	133	133			± 5	177	177	178			± 7	
	mmol/l	3.73	3.73	3.66			± 0.25	8.26	8.26	8.26			± 0.31	10.99	10.99	11.05			± 0.44	
HT HCT	%	15.9	15.9	15.8			± 2.0	36.5	36.5	36.9			± 2.5	49.7	49.7	50.4			± 3.0	
	l/l	0.159	0.159	0.158			± 0.020	0.365	0.365	0.369			± 0.025	0.497	0.497	0.504			± 0.030	
VGM MCV	µm ³ , fl	72	72	71			± 4	82	82	82			± 4	90	90	90			± 4	
TGMH MCH	pg	27.1	27.1	26.6			± 2.4	29.9	29.9	29.6			± 2.4	32.1	32.1	31.8			± 2.8	
	fmol	1.69	1.69	1.65			± 0.15	1.86	1.86	1.84			± 0.15	1.99	1.99	1.97			± 0.18	
CCMH MCHC	g/dl	37.7	37.7	37.4			± 3.0	36.4	36.4	36.0			± 3.0	35.6	35.6	35.3			± 3.0	
	g/l	377	377	374			± 30	364	364	360			± 30	356	356	353			± 30	
	mmol/l	23.42	23.42	23.25			± 1.86	22.63	22.63	22.38			± 1.86	22.13	22.13	21.93			± 1.86	
IDR RDW	%	13.7	13.7	14.0			± 3.0	13.5	13.5	13.5			± 3.0	13.0	13.0	13.0			± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	255	255	253			± 40	505	505	504			± 60	
VPM MPV	µm ³ , fl	8.3	8.3	7.8			± 2.0	7.7	7.7	7.5			± 2.0	7.7	7.7	7.5			± 2.0	
LYMPHO	#	1.4	1.4	1.3			± 0.2	2.6	2.6	2.6			± 0.5	2.9	2.9	2.9			± 1.3	
	%	64.5	64.5	64.0			± 8.0	34.0	34.0	34.0			± 6.0	14.0	14.0	14.0			± 6.0	
MONO	#	0.2	0.2	0.2			± 0.2	0.6	0.6	0.6			± 0.4	1.1	1.1	1.0			± 1.0	
	%	10.0	10.0	9.5			± 7.0	8.0	8.0	7.5			± 5.0	5.5	5.5	5.0			± 5.0	
GRANULO	#	0.5	0.5	0.6			± 0.2	4.4	4.4	4.4			± 0.5	16.4	16.4	16.5			± 1.5	
	%	25.5	25.5	26.5			± 7.0	58.0	58.0	58.5			± 6.0	80.5	80.5	81.0			± 7.0	
CRP	mg/dl	0.72	0.50	N/A			± 0.30	2.68	2.32	N/A			± 0.70	4.35	3.85	N/A			± 1.00	
	mg/l	7.2	5.0	N/A			± 3.0	26.8	23.2	N/A			± 7.0	43.5	38.5	N/A			± 10.0	

9930089-A
FRONT/RECTO
Ref: TEMP-1070 Rev.36