

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 043

CONTROL



(Exp.) 2013-07-05

(YYYY - MM - DD)

Rev 2

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE
		MICROS						MICROS						MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0			± 0.4	7.6	7.6	7.7			± 1.0	20.5	20.5	20.6			± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.35	2.33			± 0.15	4.60	4.60	4.58			± 0.18	5.72	5.72	5.72			± 0.25
HB HGB	g/dl	6.1	6.1	6.1			± 0.4	13.5	13.5	13.5			± 0.5	18.0	18.0	18.0			± 0.7
	g/l	61	61	61			± 4	135	135	135			± 5	180	180	180			± 7
	mmol/l	3.79	3.79	3.79			± 0.25	8.38	8.38	8.38			± 0.31	11.18	11.18	11.18			± 0.44
HT HCT	%	16.7	16.7	16.5			± 2.0	37.7	37.7	37.6			± 2.5	50.9	50.9	50.9			± 3.0
	l/l	0.167	0.167	0.165			± 0.020	0.377	0.377	0.376			± 0.025	0.509	0.509	0.509			± 0.030
VGM MCV	µm ³ ; fl	71	71	71			± 4	82	82	82			± 4	89	89	89			± 4
TGMH MCH	pg	26.0	26.0	26.2			± 2.4	29.3	29.3	29.5			± 2.4	31.5	31.5	31.5			± 2.8
	fmol	1.61	1.61	1.63			± 0.15	1.82	1.82	1.83			± 0.15	1.95	1.95	1.95			± 0.18
CCMH MCHC	g/dl	36.6	36.6	36.9			± 3.0	35.8	35.8	35.9			± 3.0	35.4	35.4	35.4			± 3.0
	g/l	366	366	369			± 30	358	358	359			± 30	354	354	354			± 30
	mmol/l	22.70	22.70	22.90			± 1.86	22.23	22.23	22.32			± 1.86	21.96	21.96	21.96			± 1.86
IDR RDW	%	13.6	13.6	13.8			± 3.0	13.8	13.8	13.2			± 3.0	13.0	13.0	13.0			± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	73	73	74			± 20	260	260	260			± 40	515	515	510			± 60
VPM MPV	µm ³ ; fl	8.5	8.5	8.3			± 2.0	7.8	7.8	7.5			± 2.0	7.8	7.8	7.6			± 2.0
LYMPHO	#	1.3	1.3	1.3			± 0.3	2.5	2.5	2.5			± 0.5	3.0	3.0	3.0			± 1.1
	%	64.0	64.0	62.5			± 7.0	32.3	32.3	32.0			± 6.0	14.5	14.5	14.5			± 5.0
MONO	#	0.2	0.2	0.2			± 0.2	0.6	0.6	0.6			± 0.4	1.1	1.1	1.0			± 1.0
	%	9.5	9.5	9.5			± 7.0	8.4	8.4	8.0			± 5.0	5.5	5.5	5.0			± 5.0
GRANULO	#	0.5	0.5	0.6			± 0.3	4.5	4.5	4.6			± 0.5	16.4	16.4	16.6			± 1.5
	%	26.5	26.5	28.0			± 7.0	59.3	59.3	60.0			± 6.0	80.0	80.0	80.5			± 7.0
CRP	mg/dl	1.00	0.60	N/A			± 0.30	2.80	2.10	N/A			± 0.70	4.80	4.00	N/A			± 1.00
	mg/l	10.0	6.0	N/A			± 3.0	28.0	21.0	N/A			± 7.0	48.0	40.0	N/A			± 10.0

9930074-A
FRONT/RECTO
Rev.31
Ref: TEMP-1070