

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 014

CONTROL



(Exp.) **2014-04-05**

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse													TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3			
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE			
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P				
GB WBC	$10^3/mm^3; 10^9/l$	2.1	2.1	2.0		± 0.4	7.6	7.6	7.7		± 1.0	20.2	20.2	20.5		± 2.2		
GR RBC	$10^6/mm^3; 10^{12}/l$	2.35	2.35	2.35		± 0.15	4.64	4.64	4.63		± 0.18	5.72	5.72	5.77		± 0.25		
HB HGB	g/dl	6.0	6.0	6.0		± 0.4	13.5	13.5	13.4		± 0.5	17.8	17.8	17.9		± 0.7		
	g/l	60	60	60		± 4	135	135	134		± 5	178	178	179		± 7		
	mmol/l	3.73	3.73	3.73		± 0.25	8.38	8.38	8.32		± 0.31	11.05	11.05	11.12		± 0.44		
HT HCT	%	16.2	16.2	16.2		± 2.0	37.6	37.6	37.5		± 2.5	50.3	50.3	50.8		± 3.0		
	l/l	0.162	0.162	0.162		± 0.020	0.376	0.376	0.375		± 0.025	0.503	0.503	0.508		± 0.030		
VGM MCV	$\mu m^3; fl$	69	69	69		± 4	81	81	81		± 4	88	88	88		± 4		
TGMH MCH	pg	25.5	25.5	25.5		± 2.4	29.1	29.1	28.9		± 2.4	31.1	31.1	31.0		± 2.8		
	fmol	1.59	1.59	1.59		± 0.15	1.81	1.81	1.80		± 0.15	1.93	1.93	1.93		± 0.18		
CCMH MCHC	g/dl	37.0	37.0	37.0		± 3.0	35.9	35.9	35.7		± 3.0	35.4	35.4	35.3		± 3.0		
	g/l	370	370	370		± 30	359	359	357		± 30	354	354	353		± 30		
	mmol/l	22.98	22.98	22.98		± 1.86	22.31	22.31	22.19		± 1.86	21.96	21.96	21.89		± 1.86		
IDR RDW	%	13.5	13.5	13.5		± 4.0	13.0	13.0	13.5		± 4.0	12.8	12.8	13.3		± 4.0		
PLAQ. PLTS	$10^3/mm^3; 10^9/l$	79	79	77		± 20	240	240	240		± 40	520	520	515		± 60		
VPM MPV	$\mu m^3; fl$	8.3	8.3	7.7		± 2.0	8.0	8.0	7.5		± 2.0	7.7	7.7	7.4		± 2.0		
LYMPHO	#	1.4	1.4	1.3		± 0.2	2.5	2.5	2.5		± 0.5	2.9	2.9	3.0		± 1.3		
	%	65.0	65.0	64.0		± 8.0	33.0	33.0	33.0		± 6.0	14.5	14.5	14.5		± 6.0		
MONO	#	0.2	0.2	0.2		± 0.2	0.5	0.5	0.5		± 0.4	1.3	1.3	1.1		± 1.0		
	%	10.0	10.0	10.0		± 7.0	7.0	7.0	6.5		± 5.0	6.5	6.5	5.5		± 5.0		
GRANULO	#	0.5	0.5	0.5		± 0.2	4.6	4.6	4.7		± 0.5	16.0	16.0	16.4		± 1.5		
	%	25.0	25.0	26.0		± 7.0	60.0	60.0	60.5		± 6.0	79.0	79.0	80.0		± 7.0		
CRP	mg/dl	0.92	0.60	N/A		± 0.30	2.92	2.37	N/A		± 0.70	5.01	4.49	N/A		± 1.00		
	mg/l	9.2	6.0	N/A		± 3.0	29.2	23.7	N/A		± 7.0	50.1	44.9	N/A		± 10.0		

9930074-A FRONT/RECTO Ref: TEMP-1070 Rev.33