

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 407

CONTROL



(Exp.) 2018-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse														TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3				
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE				
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0	± 0.4	7.5	7.5	7.7	± 1.0	20.1	20.1	20.5	± 2.2					
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.32	2.32	2.32	± 0.15	4.62	4.62	4.62	± 0.18	5.83	5.83	5.86	± 0.25					
HB	HGB	g/dl	6.1	6.1	6.1	± 0.4	13.3	13.3	13.4	± 0.5	17.7	17.7	17.8	± 0.7					
		g/l	61	61	61	± 4	133	133	134	± 5	177	177	178	± 7					
		mmol/l	3.79	3.79	3.79	± 0.25	8.26	8.26	8.32	± 0.31	10.99	10.99	11.05	± 0.44					
HT	HCT	%	16.2	16.2	16.5	± 2.0	37.4	37.4	37.4	± 2.5	50.1	50.1	50.4	± 3.0					
		l/l	0.162	0.162	0.165	± 0.020	0.374	0.374	0.374	± 0.025	0.501	0.501	0.504	± 0.030					
VGM	MCV	µm ³ ; fl	70	70	71	± 4	81	81	81	± 4	86	86	86	± 4					
TGMH	MCH	pg	26.3	26.3	26.3	± 2.4	28.8	28.8	29.0	± 2.4	30.4	30.4	30.4	± 2.8					
		fmol	1.63	1.63	1.63	± 0.15	1.79	1.79	1.80	± 0.15	1.89	1.89	1.89	± 0.18					
CCMH	MCHC	g/dl	37.6	37.6	37.0	± 3.0	35.5	35.5	35.8	± 3.0	35.3	35.3	35.3	± 3.0					
		g/l	376	376	370	± 30	355	355	358	± 30	353	353	353	± 30					
		mmol/l	23.33	23.33	23.00	± 1.86	22.07	22.07	22.24	± 1.86	21.92	21.92	21.93	± 1.86					
IDR	RDW	%	14.0	14.0	13.0	± 3.0	14.0	14.0	13.5	± 3.0	13.0	13.0	12.0	± 3.0					
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	79	79	81	± 20	256	256	260	± 40	509	509	510	± 60					
VPM	MPV	µm ³ ; fl	8.1	8.1	7.6	± 2.0	8.5	8.5	8.1	± 2.0	8.0	8.0	7.7	± 2.0					
LYMPHO		#	1.3	1.3	1.3	± 0.2	2.5	2.5	2.5	± 0.5	3.0	3.0	3.1	± 1.2					
		%	66.0	66.0	64.5	± 8.0	33.5	33.5	32.5	± 6.0	15.0	15.0	15.0	± 6.0					
MONO		#	0.2	0.2	0.2	± 0.2	0.5	0.5	0.6	± 0.4	0.9	0.9	0.8	± 0.8					
		%	11.0	11.0	11.0	± 7.0	7.0	7.0	7.5	± 5.0	4.5	4.5	4.0	± 4.0					
GRANULO		#	0.5	0.5	0.5	± 0.2	4.5	4.5	4.6	± 0.5	16.2	16.2	16.6	± 1.5					
		%	23.0	23.0	24.5	± 7.0	59.5	59.5	60.0	± 6.0	80.5	80.5	81.0	± 7.0					
CRP		mg/dl	0.69	0.50	N/A	± 0.30	3.14	2.65	N/A	± 0.70	4.48	3.88	N/A	± 1.00					
		mg/l	6.9	5.0	N/A	± 3.0	31.4	26.5	N/A	± 7.0	44.8	38.8	N/A	± 10.0					

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A