

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 400

CONTROL

(Exp.) **2016-10-05**
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse															
				CONTROL				1	CONTROL				2	CONTROL				3	TOLERANCES TOLERANCE
				MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS					
				CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0		± 0.4	7.3		7.3	7.5		± 1.0		19.5	19.5	19.8		± 2.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.35	2.35	2.36		± 0.15	4.60	4.60	4.63		± 0.18	5.76	5.76	5.80		± 0.25			
HB HGB	g/dl	6.2	6.2	6.1		± 0.4	13.6	13.6	13.6		± 0.5	18.1	18.1	18.1		± 0.7			
	g/l	62	62	61		± 4	136	136	136		± 5	181	181	181		± 7			
	mmol/l	3.85	3.85	3.79		± 0.25	8.45	8.45	8.45		± 0.31	11.24	11.24	11.24		± 0.44			
HT HCT	%	16.5	16.5	16.5		± 2.0	37.3	37.3	37.5		± 2.5	50.1	50.1	50.5		± 3.0			
	l/l	0.165	0.165	0.165		± 0.020	0.373	0.373	0.375		± 0.025	0.501	0.501	0.505		± 0.030			
VGM MCV	µm ³ ; fl	70	70	70		± 4	81	81	81		± 4	87	87	87		± 4			
TGMH MCH	pg	26.4	26.4	25.8		± 2.4	29.6	29.6	29.4		± 2.4	31.4	31.4	31.2		± 2.8			
	fmol	1.64	1.64	1.61		± 0.15	1.84	1.84	1.82		± 0.15	1.95	1.95	1.94		± 0.18			
CCMH MCHC	g/dl	37.7	37.7	36.9		± 3.0	36.5	36.5	36.3		± 3.0	36.1	36.1	35.9		± 3.0			
	g/l	377	377	369		± 30	365	365	363		± 30	361	361	359		± 30			
	mmol/l	23.41	23.41	22.93		± 1.86	22.67	22.67	22.52		± 1.86	22.43	22.43	22.28		± 1.86			
IDR RDW	%	13.6	13.6	13.4		± 3.0	13.2	13.2	13.3		± 3.0	12.8	12.8	13.0		± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	72		± 20	254	254	250		± 40	510	510	505		± 60			
VPM MPV	µm ³ ; fl	8.7	8.7	8.3		± 2.0	8.0	8.0	7.6		± 2.0	7.9	7.9	7.5		± 2.0			
LYMPHO	#	1.3	1.3	1.3		± 0.3	2.5	2.5	2.6		± 0.5	3.0	3.0	3.2		± 1.2			
	%	64.0	64.0	62.8		± 8.0	34.0	34.0	35.0		± 6.0	15.5	15.5	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.5		± 0.4	1.2	1.2	1.0		± 1.0			
	%	10.0	10.0	10.0		± 7.0	8.0	8.0	7.0		± 5.0	6.0	6.0	5.0		± 5.0			
GRANULO	#	0.5	0.5	0.5		± 0.3	4.2	4.2	4.4		± 0.5	15.3	15.3	15.6		± 1.4			
	%	26.0	26.0	27.2		± 8.0	58.0	58.0	58.0		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.54	0.49	N/A		± 0.30	2.65	2.65	N/A		± 0.70	4.73	4.75	N/A		± 1.00			
	mg/l	5.4	4.9	N/A		± 3.0	26.5	26.5	N/A		± 7.0	47.3	47.5	N/A		± 10.0			

Ref: TEMP-1070 Rev.36 FRONT/RECTO 9930089-A