

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MC 405**

**CONTROL**



**(Exp.) 2018-01-05**  
(YYYY - MM - DD)

Rev 1

		<b>ABX Lyse</b>																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS				MICROS		MICROS				MICROS		MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.1	2.1	2.2		± 0.4	7.6	7.6	7.6		± 1.0	20.2	20.2	20.3		± 2.2			
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.21	2.21	2.23		± 0.15	4.67	4.67	4.67		± 0.18	5.64	5.64	5.69		± 0.25			
HB HGB	g/dl	5.8	5.8	5.8		± 0.4	13.5	13.5	13.6		± 0.5	17.7	17.7	17.8		± 0.7			
	g/l	58	58	58		± 4	135	135	136		± 5	177	177	178		± 7			
HT HCT	mmol/l	3.60	3.60	3.60		± 0.25	8.38	8.38	8.45		± 0.31	10.99	10.99	11.05		± 0.44			
	%	15.5	15.5	15.6		± 2.0	37.8	37.8	37.8		± 2.5	49.6	49.6	50.6		± 3.0			
VGM MCV	l/l	0.155	0.155	0.156		± 0.020	0.378	0.378	0.378		± 0.025	0.496	0.496	0.506		± 0.030			
	µm <sup>3</sup> ; fl	70	70	70		± 4	81	81	81		± 4	88	88	89		± 4			
TGMH MCH	pg	26.2	26.2	26.0		± 2.4	28.9	28.9	29.1		± 2.4	31.4	31.4	31.3		± 2.8			
	fmol	1.63	1.63	1.62		± 0.15	1.80	1.80	1.81		± 0.15	1.95	1.95	1.94		± 0.18			
CCMH MCHC	g/dl	37.5	37.5	37.2		± 3.0	35.7	35.7	36.0		± 3.0	35.7	35.7	35.1		± 3.0			
	g/l	375	375	372		± 30	357	357	360		± 30	357	357	351		± 30			
	mmol/l	23.28	23.28	23.07		± 1.86	22.16	22.16	22.33		± 1.86	22.15	22.15	21.83		± 1.86			
IDR RDW	%	14.5	14.5	13.0		± 3.0	14.0	14.0	13.5		± 3.0	13.5	13.5	12.5		± 3.0			
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	69	69	71		± 20	255	255	257		± 40	530	530	535		± 60			
VPM MPV	µm <sup>3</sup> ; fl	8.7	8.7	8.0		± 2.0	8.1	8.1	7.7		± 2.0	8.3	8.3	7.8		± 2.0			
LYMPHO	#	1.4	1.4	1.5		± 0.2	2.7	2.7	2.7		± 0.5	3.2	3.2	3.2		± 1.3			
	%	67.5	67.5	66.0		± 8.0	35.0	35.0	35.0		± 6.0	16.0	16.0	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.5	0.5	0.5		± 0.4	1.1	1.1	1.0		± 1.0			
	%	8.5	8.5	8.5		± 7.0	7.0	7.0	6.5		± 5.0	5.5	5.5	5.0		± 5.0			
GRANULO	#	0.6	0.6	0.6		± 0.2	4.4	4.4	4.4		± 0.5	15.9	15.9	16.0		± 1.5			
	%	24.0	24.0	25.5		± 7.0	58.0	58.0	58.5		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.61	0.46	N/A		± 0.30	2.24	1.94	N/A		± 0.70	3.69	3.27	N/A		± 1.00			
	mg/l	6.1	4.6	N/A		± 3.0	22.4	19.4	N/A		± 7.0	36.9	32.7	N/A		± 10.0			

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A