

# ABX Minotrol CRP



**LOT**

**MC 406**

**CONTROL**



**(Exp.) 2018-04-05**  
(YYYY - MM - DD)

Rev 1

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS				MICROS		MICROS				MICROS		MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.1	2.1	2.1		± 0.4	7.5	7.5	7.6		± 1.0	20.1	20.1	20.4		± 2.2			
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.28	2.28	2.27		± 0.15	4.72	4.72	4.72		± 0.18	5.67	5.67	5.72		± 0.25			
HB HGB	g/dl	5.9	5.9	5.9		± 0.4	13.6	13.6	13.6		± 0.5	17.7	17.7	17.7		± 0.7			
	g/l	59	59	59		± 4	136	136	136		± 5	177	177	177		± 7			
HT HCT	mmol/l	3.66	3.66	3.66		± 0.25	8.45	8.45	8.45		± 0.31	10.99	10.99	10.99		± 0.44			
	%	16.0	16.0	16.1		± 2.0	38.2	38.2	38.2		± 2.5	49.9	49.9	50.3		± 3.0			
VGM MCV	l/l	0.160	0.160	0.161		± 0.020	0.382	0.382	0.382		± 0.025	0.499	0.499	0.503		± 0.030			
	µm <sup>3</sup> ; fl	70	70	71		± 4	81	81	81		± 4	88	88	88		± 4			
TGMH MCH	pg	25.9	25.9	26.0		± 2.4	28.8	28.8	28.8		± 2.4	31.2	31.2	30.9		± 2.8			
	fmol	1.61	1.61	1.61		± 0.15	1.79	1.79	1.79		± 0.15	1.94	1.94	1.92		± 0.18			
CCMH MCHC	g/dl	37.0	37.0	36.6		± 3.0	35.6	35.6	35.6		± 3.0	35.5	35.5	35.2		± 3.0			
	g/l	370	370	366		± 30	356	356	356		± 30	355	355	352		± 30			
	mmol/l	22.96	22.96	22.73		± 1.86	22.09	22.09	22.09		± 1.86	22.03	22.03	21.84		± 1.86			
IDR RDW	%	14.0	14.0	12.8		± 3.0	13.5	13.5	13.0		± 3.0	13.5	13.5	12.0		± 3.0			
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	76	76	80		± 20	255	255	255		± 40	495	495	505		± 60			
VPM MPV	µm <sup>3</sup> ; fl	8.3	8.3	7.8		± 2.0	8.5	8.5	8.1		± 2.0	7.8	7.8	7.4		± 2.0			
LYMPHO	#	1.3	1.3	1.3		± 0.2	2.7	2.7	2.7		± 0.5	3.3	3.3	3.3		± 1.3			
	%	64.0	64.0	63.5		± 8.0	36.0	36.0	35.0		± 6.0	16.5	16.5	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.5	0.5	0.5		± 0.5	1.0	1.0	1.0		± 1.0			
	%	8.5	8.5	9.0		± 7.0	6.0	6.0	6.0		± 5.0	5.0	5.0	5.0		± 5.0			
GRANULO	#	0.6	0.6	0.6		± 0.2	4.4	4.4	4.5		± 0.5	15.8	15.8	16.1		± 1.5			
	%	27.5	27.5	27.5		± 7.0	58.0	58.0	59.0		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.80	0.57	N/A		± 0.30	2.55	2.14	N/A		± 0.70	4.81	4.28	N/A		± 1.00			
	mg/l	8.0	5.7	N/A		± 3.0	25.5	21.4	N/A		± 7.0	48.1	42.8	N/A		± 10.0			

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A