

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 409

CONTROL



(Exp.) **2019-01-05**
(YYYY - MM - DD)

Rev 1

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS				MICROS		MICROS				MICROS		MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1	2.1		± 0.4	7.5	7.5	7.6		± 1.0	19.9	19.9	20.3		± 2.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.29	2.29	2.32		± 0.15	4.63	4.63	4.64		± 0.18	5.85	5.85	5.87		± 0.25			
HB HGB	g/dl	6.1	6.1	6.2		± 0.4	13.4	13.4	13.5		± 0.5	18.1	18.1	18.2		± 0.7			
	g/l	61	61	62		± 4	134	134	135		± 5	181	181	182		± 7			
HT HCT	mmol/l	3.79	3.79	3.85		± 0.25	8.32	8.32	8.38		± 0.31	11.24	11.24	11.30		± 0.44			
	%	16.3	16.3	16.7		± 2.0	37.5	37.5	37.6		± 2.5	51.5	51.5	51.7		± 3.0			
VGM MCV	l/l	0.163	0.163	0.167		± 0.020	0.375	0.375	0.376		± 0.025	0.515	0.515	0.517		± 0.030			
	µm ³ ; fl	71	71	72		± 4	81	81	81		± 4	88	88	88		± 4			
TGMH MCH	pg	26.6	26.6	26.7		± 2.4	28.9	28.9	29.1		± 2.4	30.9	30.9	31.0		± 2.8			
	fmol	1.65	1.65	1.66		± 0.15	1.80	1.80	1.81		± 0.15	1.92	1.92	1.93		± 0.18			
CCMH MCHC	g/dl	37.5	37.5	37.1		± 3.0	35.7	35.7	35.9		± 3.0	35.2	35.2	35.2		± 3.0			
	g/l	375	375	371		± 30	357	357	359		± 30	352	352	352		± 30			
	mmol/l	23.30	23.30	23.05		± 1.86	22.19	22.19	22.31		± 1.86	21.83	21.83	21.88		± 1.86			
IDR RDW	%	13.5	13.5	13.0		± 3.0	13.5	13.5	13.0		± 3.0	13.0	13.0	12.0		± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	76	76	80		± 20	255	255	257		± 40	500	500	505		± 60			
VPM MPV	µm ³ ; fl	8.3	8.3	7.5		± 2.0	8.1	8.1	7.3		± 2.0	7.9	7.9	7.2		± 2.0			
LYMPHO	#	1.4	1.4	1.3		± 0.2	2.6	2.6	2.5		± 0.5	3.2	3.2	3.2		± 1.2			
	%	65.0	65.0	63.0		± 8.0	34.0	34.0	33.5		± 6.0	16.0	16.0	16.0		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6		± 0.4	1.1	1.1	1.0		± 1.0			
	%	9.0	9.0	9.5		± 7.0	8.0	8.0	7.5		± 5.0	5.5	5.5	5.0		± 5.0			
GRANULO	#	0.5	0.5	0.6		± 0.2	4.4	4.4	4.5		± 0.6	15.6	15.6	16.0		± 1.5			
	%	26.0	26.0	27.5		± 7.0	58.0	58.0	59.0		± 6.0	78.5	78.5	79.0		± 7.0			
CRP	mg/dl	0.65	0.55	N/A		± 0.30	2.80	2.70	N/A		± 0.70	5.10	4.89	N/A		± 1.00			
	mg/l	6.5	5.5	N/A		± 3.0	28.0	27.0	N/A		± 7.0	51.0	48.9	N/A		± 10.0			

Ref: TEMP-1070 Rev.37 FRONT/RECTO 9930089-A