

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 103

CONTROL



(Exp.) 2014-01-05

(YYYY - MM - DD)

Rev 1

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE	CONTROL					TOLERANCES TOLERANCE
		MICROS						MICROS						MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0			± 0.4	7.4	7.4	7.6			± 1.0	19.8	19.8	20.3			± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.39	2.39	2.39			± 0.15	4.57	4.57	4.59			± 0.18	5.53	5.53	5.59			± 0.25
HB HGB	g/dl	6.1	6.1	6.1			± 0.4	13.4	13.4	13.4			± 0.5	17.5	17.5	17.5			± 0.7
	g/l	61	61	61			± 4	134	134	134			± 5	175	175	175			± 7
	mmol/l	3.79	3.79	3.79			± 0.25	8.32	8.32	8.32			± 0.31	10.87	10.87	10.87			± 0.44
HT HCT	%	16.7	16.7	16.7			± 2.0	37.0	37.0	37.2			± 2.5	49.2	49.2	49.7			± 3.0
	l/l	0.167	0.167	0.167			± 0.020	0.370	0.370	0.372			± 0.025	0.492	0.492	0.497			± 0.030
VGM MCV	µm ³ ; fl	70	70	70			± 4	81	81	81			± 4	89	89	89			± 4
TGMH MCH	pg	25.5	25.5	25.5			± 2.4	29.3	29.3	29.2			± 2.4	31.6	31.6	31.3			± 2.8
	fmol	1.58	1.58	1.58			± 0.15	1.82	1.82	1.81			± 0.15	1.97	1.97	1.94			± 0.18
CCMH MCHC	g/dl	36.5	36.5	36.5			± 3.0	36.2	36.2	36.0			± 3.0	35.6	35.6	35.2			± 3.0
	g/l	365	365	365			± 30	362	362	360			± 30	356	356	352			± 30
	mmol/l	22.64	22.64	22.64			± 1.86	22.48	22.48	22.38			± 1.86	22.08	22.08	21.84			± 1.86
IDR RDW	%	14.0	14.0	14.0			± 3.0	13.5	13.5	13.0			± 3.0	12.5	12.5	12.5			± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	75	75			± 20	240	240	246			± 40	490	490	495			± 60
VPM MPV	µm ³ ; fl	8.7	8.7	8.1			± 2.0	7.7	7.7	7.3			± 2.0	7.6	7.6	7.3			± 2.0
LYMPHO	#	1.3	1.3	1.3			± 0.2	2.4	2.4	2.5			± 0.5	3.0	3.0	3.0			± 1.2
	%	64.0	64.0	62.5			± 8.0	32.5	32.5	32.5			± 6.0	15.0	15.0	15.0			± 6.0
MONO	#	0.2	0.2	0.2			± 0.2	0.7	0.7	0.6			± 0.4	1.4	1.4	1.2			± 1.0
	%	10.0	10.0	10.5			± 7.0	9.5	9.5	8.5			± 5.0	7.0	7.0	6.0			± 5.0
GRANULO	#	0.5	0.5	0.5			± 0.2	4.3	4.3	4.5			± 0.5	15.4	15.4	16.0			± 1.5
	%	26.0	26.0	27.0			± 7.0	58.0	58.0	59.0			± 6.0	78.0	78.0	79.0			± 7.0
CRP	mg/dl	0.84	0.50	N/A			± 0.30	2.52	1.93	N/A			± 0.70	4.44	3.71	N/A			± 1.00
	mg/l	8.4	5.0	N/A			± 3.0	25.2	19.3	N/A			± 7.0	44.4	37.1	N/A			± 10.0

9930074-A FRONT/RECTO Ref: TEMP-1070 Rev.31