

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 044

CONTROL



(Exp.) 2014-07-05

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse													TOLERANCES TOLERANCE		
				CONTROL			1	CONTROL			2		CONTROL			3			
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE		MICROS			TOLERANCES TOLERANCE			
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P			CRP	CRP 200	60 8P 60 18P				
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.0		± 0.4	7.4	7.4	7.5			± 1.0	20.1	20.1	20.5			± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.37	2.37	2.39		± 0.15	4.60	4.60	4.64			± 0.18	5.66	5.66	5.72			± 0.25
HB	HGB	g/dl	6.1	6.1	6.0		± 0.4	13.5	13.5	13.5			± 0.5	17.9	17.9	17.9			± 0.7
		g/l	61	61	60		± 4	135	135	135			± 5	179	179	179			± 7
		mmol/l	3.79	3.79	3.73		± 0.25	8.38	8.38	8.38			± 0.31	11.12	11.12	11.12			± 0.44
HT	HCT	%	16.4	16.4	16.5		± 2.0	37.3	37.3	37.6			± 2.5	50.4	50.4	50.9			± 3.0
		l/l	0.164	0.164	0.165		± 0.020	0.373	0.373	0.376			± 0.025	0.504	0.504	0.509			± 0.030
VGM	MCV	µm ³ ; fl	69	69	69		± 4	81	81	81			± 4	89	89	89			± 4
TGMH	MCH	pg	25.7	25.7	25.1		± 2.4	29.3	29.3	29.1			± 2.4	31.6	31.6	31.3			± 2.8
		fmol	1.60	1.60	1.56		± 0.15	1.82	1.82	1.81			± 0.15	1.96	1.96	1.94			± 0.18
CCMH	MCHC	g/dl	37.3	37.3	36.4		± 3.0	36.2	36.2	35.9			± 3.0	35.5	35.5	35.2			± 3.0
		g/l	373	373	364		± 30	362	362	359			± 30	355	355	352			± 30
		mmol/l	23.16	23.16	22.59		± 1.86	22.50	22.50	22.31			± 1.86	22.07	22.07	21.84			± 1.86
IDR	RDW	%	13.5	13.5	13.5		± 3.0	13.4	13.4	13.6			± 3.0	13.2	13.2	13.5			± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	79	79	79		± 20	255	255	255			± 40	500	500	505			± 60
VPM	MPV	µm ³ ; fl	8.8	8.8	8.4		± 2.0	8.1	8.1	7.7			± 2.0	7.8	7.8	7.5			± 2.0
LYMPHO		#	1.3	1.3	1.3		± 0.3	2.5	2.5	2.6			± 0.5	3.0	3.0	3.1			± 1.3
		%	65.5	65.5	64.5		± 8.0	34.0	34.0	34.0			± 6.0	15.0	15.0	15.0			± 6.0
MONO		#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6			± 0.4	1.4	1.4	1.1			± 1.0
		%	9.5	9.5	9.0		± 7.0	8.5	8.5	7.5			± 5.0	7.0	7.0	5.5			± 5.0
GRANULO		#	0.5	0.5	0.5		± 0.2	4.3	4.3	4.4			± 0.5	15.7	15.7	16.3			± 1.5
		%	25.0	25.0	26.5		± 7.0	57.5	57.5	58.5			± 6.0	78.0	78.0	79.5			± 7.0
CRP		mg/dl	0.90	0.50	N/A		± 0.30	2.90	2.50	N/A			± 0.70	4.60	4.10	N/A			± 1.00
		mg/l	9.0	5.0	N/A		± 3.0	29.0	25.0	N/A			± 7.0	46.0	41.0	N/A			± 10.0

9930074-A FRONT/RECTO Ref: TEMP-1070 Rev.34