

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 074

CONTROL



(Exp.) 2014-10-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse														TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3				
				MICROS				MICROS				MICROS							
				CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P		CRP	CRP 200	60 8P 60 18P					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.1		± 0.4	7.4	7.4	7.6		± 1.0	19.7	19.7	20.1		± 2.2		
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.34	2.34	2.36		± 0.15	4.59	4.59	4.64		± 0.18	5.47	5.47	5.53		± 0.25		
HB	HGB	g/dl	6.0	6.0	6.0		± 0.4	13.4	13.4	13.4		± 0.5	17.7	17.7	17.7		± 0.7		
		g/l	60	60	60		± 4	134	134	134		± 5	177	177	177		± 7		
		mmol/l	3.73	3.73	3.73		± 0.25	8.32	8.32	8.32		± 0.31	10.99	10.99	10.99		± 0.44		
HT	HCT	%	16.4	16.4	16.3		± 2.0	37.2	37.2	37.6		± 2.5	49.2	49.2	49.8		± 3.0		
		l/l	0.164	0.164	0.163		± 0.020	0.372	0.372	0.376		± 0.025	0.492	0.492	0.498		± 0.030		
VGM	MCV	µm ³ ; fl	70	70	69		± 4	81	81	81		± 4	90	90	90		± 4		
TGMH	MCH	pg	25.6	25.6	25.4		± 2.4	29.2	29.2	28.9		± 2.4	32.4	32.4	32.0		± 2.8		
		fmol	1.59	1.59	1.58		± 0.15	1.81	1.81	1.79		± 0.15	2.01	2.01	1.99		± 0.18		
CCMH	MCHC	g/dl	36.6	36.6	36.8		± 3.0	36.0	36.0	35.7		± 3.0	36.0	36.0	35.6		± 3.0		
		g/l	366	366	368		± 30	360	360	357		± 30	360	360	356		± 30		
		mmol/l	22.75	22.75	22.88		± 1.86	22.38	22.38	22.14		± 1.86	22.33	22.33	22.08		± 1.86		
IDR	RDW	%	14.0	14.0	14.0		± 3.0	13.7	13.7	13.7		± 3.0	13.0	13.0	13.2		± 3.0		
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	80	80	80		± 20	257	257	258		± 40	495	495	500		± 60		
VPM	MPV	µm ³ ; fl	8.6	8.6	8.0		± 2.0	7.7	7.7	7.4		± 2.0	7.6	7.6	7.3		± 2.0		
LYMPHO		#	1.3	1.3	1.3		± 0.2	2.5	2.5	2.5		± 0.5	3.1	3.1	3.1		± 1.1		
		%	63.5	63.5	63.0		± 8.0	33.5	33.5	33.5		± 6.0	15.5	15.5	15.5		± 5.0		
MONO		#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6		± 0.4	1.0	1.0	1.0		± 1.0		
		%	8.5	8.5	8.5		± 7.0	8.0	8.0	7.5		± 5.0	5.0	5.0	5.0		± 5.0		
GRANULO		#	0.6	0.6	0.6		± 0.2	4.3	4.3	4.5		± 0.5	15.7	15.7	16.0		± 1.5		
		%	28.0	28.0	28.5		± 7.0	58.5	58.5	59.0		± 6.0	79.5	79.5	79.5		± 7.0		
CRP		mg/dl	0.81	0.46	N/A		± 0.30	2.84	2.32	N/A		± 0.70	4.46	3.85	N/A		± 1.00		
		mg/l	8.1	4.6	N/A		± 3.0	28.4	23.2	N/A		± 7.0	44.6	38.5	N/A		± 10.0		

9930089-A
FRONT/RECTO
Ref: TEMP-1070 Rev.35