

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 045

CONTROL



(Exp.) 2015-07-05
(YYYY - MM - DD)

Rev 2

		ABX Lyse																	
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
		MICROS						MICROS						MICROS					
		CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	1.9	1.9	1.9		± 0.4	7.6	7.6	7.7		± 1.0	20.2	20.2	20.5		± 2.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.52	2.52	2.56		± 0.15	4.61	4.61	4.65		± 0.18	5.77	5.77	5.86		± 0.25			
HB HGB	g/dl	6.2	6.2	6.2		± 0.4	13.4	13.4	13.5		± 0.5	18.0	18.0	18.1		± 0.7			
	g/l	62	62	62		± 4	134	134	135		± 5	180	180	181		± 7			
HT HCT	mmol/l	3.85	3.85	3.85		± 0.25	8.32	8.32	8.38		± 0.31	11.18	11.18	11.24		± 0.44			
	%	17.4	17.4	17.4		± 2.0	37.3	37.3	37.7		± 2.5	51.4	51.4	51.6		± 3.0			
VGM MCV	l/l	0.174	0.174	0.174		± 0.020	0.373	0.373	0.377		± 0.025	0.514	0.514	0.516		± 0.030			
	µm ³ ; fl	69	69	68		± 4	81	81	81		± 4	89	89	88		± 4			
TGMH MCH	pg	24.6	24.6	24.2		± 2.4	29.1	29.1	29.0		± 2.4	31.2	31.2	30.9		± 2.8			
	fmol	1.53	1.53	1.50		± 0.15	1.81	1.81	1.80		± 0.15	1.94	1.94	1.92		± 0.18			
CCMH MCHC	g/dl	35.7	35.7	35.6		± 3.0	35.9	35.9	35.8		± 3.0	35.1	35.1	35.1		± 3.0			
	g/l	357	357	356		± 30	359	359	358		± 30	351	351	351		± 30			
	mmol/l	22.14	22.14	22.12		± 1.86	22.28	22.28	22.26		± 1.86	21.77	21.77	21.80		± 1.86			
IDR RDW	%	14.0	14.0	14.0		± 3.0	13.5	13.5	13.5		± 3.0	13.0	13.0	13.0		± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	70	70	75		± 20	255	255	257		± 40	505	505	508		± 60			
VPM MPV	µm ³ ; fl	9.0	9.0	8.4		± 2.0	7.7	7.7	7.4		± 2.0	7.6	7.6	7.4		± 2.0			
LYMPHO	#	1.3	1.3	1.2		± 0.3	2.6	2.6	2.6		± 0.5	3.2	3.2	3.2		± 1.3			
	%	66.0	66.0	65.0		± 8.0	34.0	34.0	33.2		± 6.0	16.0	16.0	15.5		± 6.0			
MONO	#	0.2	0.2	0.2		± 0.2	0.6	0.6	0.6		± 0.5	1.3	1.3	1.1		± 1.0			
	%	10.0	10.0	10.0		± 7.0	8.0	8.0	8.0		± 5.0	6.5	6.5	5.5		± 5.0			
GRANULO	#	0.5	0.5	0.5		± 0.2	4.4	4.4	4.5		± 0.6	15.7	15.7	16.2		± 1.5			
	%	24.0	24.0	25.0		± 7.0	58.0	58.0	58.8		± 6.0	77.5	77.5	79.0		± 7.0			
CRP	mg/dl	0.69	0.46	N/A		± 0.30	2.79	2.47	N/A		± 0.70	4.40	4.00	N/A		± 1.00			
	mg/l	6.9	4.6	N/A		± 3.0	27.9	24.7	N/A		± 7.0	44.0	40.0	N/A		± 10.0			

Ref: TEMP-1070 Rev.36 FRONT/RECTO 9930089-A