

ABX Minotrol CRP

LOT
MC 105
Rev 2
CONTROL

(Exp.) 2016-01-05

(YYYY - MM - DD)

		REV 2	ABX Lyse																(TTTT - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS		CONTROL				1	TOLERANCES TOLERANCE	CONTROL				2	TOLERANCES TOLERANCE	CONTROL				3	TOLERANCES TOLERANCE
			MICROS						MICROS						MICROS					
			CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P				CRP	CRP 200	60 8P 60 18P			
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1	2.1			± 0.4	7.9	7.9	7.9			± 1.0	19.9	19.9	19.9			± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.40	2.40	2.43			± 0.15	4.53	4.53	4.55			± 0.18	5.62	5.62	5.65			± 0.25
HB	HGB	g/dl	6.1	6.1	6.1			± 0.4	13.3	13.3	13.4			± 0.5	17.5	17.5	17.8			± 0.7
		g/l	61	61	61			± 4	133	133	134			± 5	175	175	178			± 7
		mmol/l	3.79	3.79	3.79			± 0.25	8.26	8.26	8.32			± 0.31	10.87	10.87	11.05			± 0.44
HT	HCT	%	16.8	16.8	16.8			± 2.0	37.1	37.1	37.3			± 2.5	50.0	50.0	50.3			± 3.0
		l/l	0.168	0.168	0.168			± 0.020	0.371	0.371	0.373			± 0.025	0.500	0.500	0.503			± 0.030
VGM	MCV	µm ³ ; fl	70	70	69			± 4	82	82	82			± 4	89	89	89			± 4
TGMH	MCH	pg	25.4	25.4	25.1			± 2.4	29.4	29.4	29.5			± 2.4	31.1	31.1	31.5			± 2.8
		fmol	1.58	1.58	1.56			± 0.15	1.82	1.82	1.83			± 0.15	1.93	1.93	1.96			± 0.18
CCMH	MCHC	g/dl	36.3	36.3	36.4			± 3.0	35.8	35.8	35.9			± 3.0	35.0	35.0	35.4			± 3.0
		g/l	363	363	364			± 30	358	358	359			± 30	350	350	354			± 30
		mmol/l	22.55	22.55	22.59			± 1.86	22.23	22.23	22.30			± 1.86	21.73	21.73	21.98			± 1.86
IDR	RDW	%	13.8	13.8	13.8			± 3.0	13.5	13.5	13.3			± 3.0	12.7	12.7	13.0			± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	66	66	71			± 20	248	248	250			± 40	480	480	480			± 60
VPM	MPV	µm ³ ; fl	9.2	9.2	8.5			± 2.0	8.0	8.0	7.7			± 2.0	7.8	7.8	7.3			± 2.0
LYMPHO	#	1.3	1.3	1.3			± 0.3	2.8	2.8	2.7			± 0.5	3.2	3.2	3.1			± 1.3	
	%	64.0	64.0	62.5			± 8.0	35.0	35.0	34.0			± 6.0	16.0	16.0	15.5			± 6.0	
MONO	#	0.2	0.2	0.2			± 0.2	0.7	0.7	0.7			± 0.5	1.2	1.2	1.1			± 1.0	
	%	11.5	11.5	11.0			± 7.0	9.0	9.0	8.5			± 5.0	6.0	6.0	5.5			± 5.0	
GRANULO	#	0.5	0.5	0.6			± 0.2	4.4	4.4	4.5			± 0.6	15.5	15.5	15.7			± 1.5	
	%	24.5	24.5	26.5			± 7.0	56.0	56.0	57.5			± 6.0	78.0	78.0	79.0			± 7.0	
CRP	mg/dl	0.78	0.62	N/A			± 0.30	2.40	2.21	N/A			± 0.70	4.13	4.01	N/A			± 1.00	
	mg/l	7.8	6.2	N/A			± 3.0	24.0	22.1	N/A			± 7.0	41.3	40.1	N/A			± 10.0	

Ref: TEMP-1070 Rev.36 FRONT/RECTO 9930089-A