

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 014

CONTROL



(Exp.) **2014-04-05**

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		Minilysebio																													
		CONTROL					1					CONTROL					2					CONTROL					3				
		MICROS					TOLERANCES					MICROS					TOLERANCES					MICROS					TOLERANCES				
		EMI CRP	CRP 200				TOLERANCE	EMI CRP	CRP 200				TOLERANCE	EMI CRP	CRP 200				TOLERANCE	EMI CRP	CRP 200				TOLERANCE						
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1				± 0.4	7.8	7.6					± 1.0	20.4	19.8							± 2.2							
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.37	2.35				± 0.15	4.60	4.64					± 0.18	5.62	5.72							± 0.25							
HB	HGB	g/dl	6.0	5.7				± 0.4	13.3	13.2					± 0.5	17.7	17.8							± 0.7							
		g/l	60	57				± 4	133	132					± 5	177	178							± 7							
		mmol/l	3.73	3.54				± 0.25	8.26	8.22					± 0.31	10.99	11.05							± 0.44							
HT	HCT	%	16.4	16.2				± 2.0	37.3	37.6					± 2.5	49.5	50.3							± 3.0							
		l/l	0.164	0.162				± 0.020	0.373	0.376					± 0.025	0.495	0.503							± 0.030							
VGM	MCV	µm ³ ; fl	69	69				± 4	81	81					± 4	88	88							± 4							
TGMH	MCH	pg	25.3	24.3				± 2.4	28.9	28.5					± 2.4	31.5	31.1							± 2.8							
		fmol	1.57	1.51				± 0.15	1.80	1.77					± 0.15	1.96	1.93							± 0.18							
CCMH	MCHC	g/dl	36.7	35.2				± 3.0	35.7	35.2					± 3.0	35.8	35.4							± 3.0							
		g/l	367	352				± 30	357	352					± 30	358	354							± 30							
		mmol/l	22.78	21.83				± 1.86	22.17	21.86					± 1.86	22.23	21.96							± 1.86							
IDR	RDW	%	13.2	13.5				± 4.0	12.8	13.0					± 4.0	11.9	12.8							± 4.0							
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	79				± 20	235	240					± 40	505	520							± 60							
VPM	MPV	µm ³ ; fl	8.3	8.3				± 2.0	8.5	8.0					± 2.0	8.2	7.7							± 2.0							
LYMPHO		#	1.3	1.4				± 0.2	2.6	2.6					± 0.5	3.4	3.1							± 1.3							
		%	63.0	67.6				± 8.0	33.5	34.7					± 6.0	16.5	15.5							± 6.0							
MONO		#	0.2	0.2				± 0.2	0.5	0.5					± 0.4	1.1	1.3							± 1.0							
		%	9.0	8.0				± 7.0	6.0	6.7					± 5.0	5.5	6.5							± 5.0							
GRANULO		#	0.6	0.5				± 0.2	4.7	4.5					± 0.5	15.9	15.5							± 1.5							
		%	28.0	24.5				± 7.0	60.5	58.8					± 6.0	78.0	78.2							± 7.0							
CRP		mg/dl	0.51	0.60				± 0.30	2.28	2.37					± 0.70	4.48	4.49							± 1.00							
		mg/l	5.1	6.0				± 3.0	22.8	23.7					± 7.0	44.8	44.9							± 10.0							

9930087-B
Ref: TEMP-1070 Rev.33