

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT**

**MC 104**

**CONTROL**



**(Exp.) 2015-01-05**

(YYYY - MM - DD)

Rev 1

		<b>Minilysebio</b>															
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	
		MICROS					MICROS					MICROS					
		CRP 200					CRP 200					CRP 200					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>2.1</b>				± 0.4	<b>7.6</b>					± 1.0	<b>20.0</b>				± 2.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>2.21</b>				± 0.15	<b>4.45</b>					± 0.18	<b>5.52</b>				± 0.25
HB HGB	g/dl	<b>5.7</b>				± 0.4	<b>13.0</b>					± 0.5	<b>17.7</b>				± 0.7
	g/l	<b>57</b>				± 4	<b>130</b>					± 5	<b>177</b>				± 7
	mmol/l	<b>3.54</b>				± 0.25	<b>8.09</b>					± 0.31	<b>10.99</b>				± 0.44
HT HCT	%	<b>15.9</b>				± 2.0	<b>36.5</b>					± 2.5	<b>49.7</b>				± 3.0
	l/l	<b>0.159</b>				± 0.020	<b>0.365</b>					± 0.025	<b>0.497</b>				± 0.030
VGM MCV	µm <sup>3</sup> , fl	<b>72</b>				± 4	<b>82</b>					± 4	<b>90</b>				± 4
TGMH MCH	pg	<b>25.8</b>				± 2.4	<b>29.3</b>					± 2.4	<b>32.1</b>				± 2.8
	fmol	<b>1.60</b>				± 0.15	<b>1.82</b>					± 0.15	<b>1.99</b>				± 0.18
CCMH MCHC	g/dl	<b>35.8</b>				± 3.0	<b>35.7</b>					± 3.0	<b>35.6</b>				± 3.0
	g/l	<b>358</b>				± 30	<b>357</b>					± 30	<b>356</b>				± 30
	mmol/l	<b>22.25</b>				± 1.86	<b>22.18</b>					± 1.86	<b>22.13</b>				± 1.86
IDR RDW	%	<b>13.7</b>				± 3.0	<b>13.5</b>					± 3.0	<b>13.0</b>				± 3.0
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>75</b>				± 20	<b>255</b>					± 40	<b>505</b>				± 60
VPM MPV	µm <sup>3</sup> , fl	<b>8.3</b>				± 2.0	<b>7.7</b>					± 2.0	<b>7.7</b>				± 2.0
LYMPHO	#	<b>1.4</b>				± 0.2	<b>2.7</b>					± 0.5	<b>3.0</b>				± 1.3
	%	<b>67.1</b>				± 8.0	<b>35.7</b>					± 6.0	<b>15.0</b>				± 6.0
MONO	#	<b>0.2</b>				± 0.2	<b>0.6</b>					± 0.4	<b>1.1</b>				± 1.0
	%	<b>8.0</b>				± 7.0	<b>7.6</b>					± 5.0	<b>5.5</b>				± 5.0
GRANULO	#	<b>0.5</b>				± 0.2	<b>4.3</b>					± 0.5	<b>15.9</b>				± 1.5
	%	<b>25.0</b>				± 7.0	<b>56.8</b>					± 6.0	<b>79.7</b>				± 7.0
CRP	mg/dl	<b>0.50</b>				± 0.30	<b>2.32</b>					± 0.70	<b>3.85</b>				± 1.00
	mg/l	<b>5.0</b>				± 3.0	<b>23.2</b>					± 7.0	<b>38.5</b>				± 10.0

9930093-A  
Ref: TEMP-1070 Rev.36