

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 015

CONTROL



(Exp.) **2015-04-05**

(YYYY - MM - DD)

Rev 1

		Minilysebio														
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE
		MICROS					MICROS					MICROS				
		CRP 200					CRP 200					CRP 200				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1				± 0.4	7.6				± 1.0	19.9				± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.34				± 0.15	4.60				± 0.18	5.60				± 0.25
HB HGB	g/dl	5.8				± 0.4	13.4				± 0.5	17.9				± 0.7
	g/l	58				± 4	134				± 5	179				± 7
	mmol/l	3.60				± 0.25	8.34				± 0.31	11.12				± 0.44
HT HCT	%	16.6				± 2.0	37.7				± 2.5	49.8				± 3.0
	l/l	0.166				± 0.020	0.377				± 0.025	0.498				± 0.030
VGM MCV	µm ³ ; fl	71				± 4	82				± 4	89				± 4
TGMH MCH	pg	24.8				± 2.4	29.2				± 2.4	32.0				± 2.8
	fmol	1.54				± 0.15	1.81				± 0.15	1.98				± 0.18
CCMH MCHC	g/dl	34.9				± 3.0	35.6				± 3.0	35.9				± 3.0
	g/l	349				± 30	356				± 30	359				± 30
	mmol/l	21.66				± 1.86	22.10				± 1.86	22.30				± 1.86
IDR RDW	%	14.0				± 3.0	13.5				± 3.0	13.0				± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	74				± 20	250				± 40	500				± 60
VPM MPV	µm ³ ; fl	8.5				± 2.0	8.0				± 2.0	7.7				± 2.0
LYMPHO	#	1.5				± 0.2	2.7				± 0.5	3.1				± 1.3
	%	69.7				± 8.0	35.7				± 6.0	15.5				± 6.0
MONO	#	0.2				± 0.2	0.5				± 0.4	1.2				± 1.0
	%	7.2				± 7.0	7.1				± 5.0	6.0				± 5.0
GRANULO	#	0.5				± 0.2	4.4				± 0.5	15.7				± 1.5
	%	23.5				± 7.0	57.3				± 6.0	78.7				± 7.0
CRP	mg/dl	0.41				± 0.30	2.12				± 0.70	4.21				± 1.00
	mg/l	4.1				± 3.0	21.2				± 7.0	42.1				± 10.0

9930093-A
Ref: TEMP-1070 Rev.36