

ABX Minotrol CRP



LOT

MC 075

Rev 2

CONTROL



(Exp.) **2015-10-05**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio															
		CONTROL				1	CONTROL				2	CONTROL				3	
		MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	
			CRP 200					CRP 200					CRP 200				
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l		2,0			± 0,4		7,6				± 1,0		19,8			± 2,2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2,41			± 0,15		4,60				± 0,18		5,68			± 0,25
HB HGB	g/dl		5,8			± 0,4		13,2				± 0,5		17,8			± 0,7
	g/l		58			± 4		132				± 5		178			± 7
	mmol/l		3,60			± 0,25		8,22				± 0,31		11,05			± 0,44
HT HCT	%		16,6			± 2,0		37,3				± 2,5		50,0			± 3,0
	l/l		0,166			± 0,020		0,373				± 0,025		0,500			± 0,030
VGM MCV	µm ³ ; fl		69			± 4		81				± 4		88			± 4
TGMH MCH	pg		24,0			± 2,4		28,8				± 2,4		31,3			± 2,8
	fmol		1,49			± 0,15		1,79				± 0,15		1,95			± 0,18
CCMH MCHC	g/dl		34,8			± 3,0		35,5				± 3,0		35,6			± 3,0
	g/l		348			± 30		355				± 30		356			± 30
	mmol/l		21,64			± 1,86		22,05				± 1,86		22,11			± 1,86
IDR RDW	%		13,8			± 3,0		13,5				± 3,0		13,2			± 3,0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		72			± 20		246				± 40		510			± 60
VPM MPV	µm ³ ; fl		9,0			± 2,0		8,3				± 2,0		8,0			± 2,0
LYMPHO	#		1,3			± 0,3		2,7				± 0,5		3,2			± 1,3
	%		67,1			± 8,0		35,7				± 6,0		16,1			± 6,0
MONO	#		0,2			± 0,2		0,6				± 0,5		1,2			± 1,0
	%		9,2			± 7,0		8,1				± 5,0		6,0			± 5,0
GRANULO	#		0,5			± 0,2		4,3				± 0,6		15,5			± 1,5
	%		23,5			± 7,0		56,4				± 6,0		78,2			± 7,0
CRP	mg/dl		0,59			± 0,30		2,09				± 0,70		4,23			± 1,00
	mg/l		5,9			± 3,0		20,9				± 7,0		42,3			± 10,0

9930093-A
Ref: TEMP-1070 Rev.36