

ABX Minotrol CRP



LOT

MC 402

Rev 2

CONTROL



(Exp.) **2017-04-05**

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio																			
				CONTROL				1		CONTROL				2		CONTROL				3		TOLERANCES TOLERANCE	
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA					
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0		2.0		± 0.4	7.8		8.1		± 1.0	20.6		21.5		± 2.2						
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.43		2.42		± 0.15	4.67		4.70		± 0.18	5.51		5.61		± 0.25						
HB	HGB	g/dl	6.3		6.0		± 0.4	13.6		13.7		± 0.5	17.5		17.9		± 0.7						
		g/l	63		60		± 4	136		137		± 5	175		179		± 7						
		mmol/l	3.91		3.73		± 0.25	8.45		8.51		± 0.31	10.87		11.12		± 0.44						
HT	HCT	%	17.0		16.9		± 2.0	37.4		37.1		± 2.5	48.5		48.2		± 3.0						
		l/l	0.170		0.169		± 0.020	0.374		0.371		± 0.025	0.485		0.482		± 0.030						
VGM	MCV	µm ³ ; fl	70		70.0		± 4	80		79.0		± 4	88		86.0		± 4						
TGMH	MCH	pg	25.9		24.8		± 2.4	29.1		29.1		± 2.4	31.8		31.9		± 2.8						
		fmol	1.61		1.54		± 0.15	1.81		1.81		± 0.15	1.97		1.98		± 0.18						
CCMH	MCHC	g/dl	37.0		35.4		± 3.0	36.4		36.9		± 3.0	36.1		37.1		± 3.0						
		g/l	370		354		± 30	364		369		± 30	361		371		± 30						
		mmol/l	23.00		22.00		± 1.86	22.61		22.91		± 1.86	22.41		23.04		± 1.86						
IDR	RDW	%	13.0		13.4		± 3.0	12.5		12.5		± 3.0	12.0		12.0		± 3.0						
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	70		65		± 20	250		245		± 40	505		495		± 60						
VPM	MPV	µm ³ ; fl	8.0		9.6		± 2.0	8.6		9.4		± 2.0	8.7		9.3		± 2.0						
LYMPHO	#		1.3		N/A		± 0.3	2.8		N/A		± 0.6	3.8		N/A		± 1.3						
	%		63.0		N/A		± 8.0	36.0		N/A		± 6.0	18.5		N/A		± 6.0						
MONO	#		0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0						
	%		9.0		N/A		± 7.0	6.0		N/A		± 5.0	5.0		N/A		± 5.0						
GRANULO	#		0.6		N/A		± 0.3	4.5		N/A		± 0.6	15.8		N/A		± 1.5						
	%		28.0		N/A		± 8.0	58.0		N/A		± 6.0	76.5		N/A		± 7.0						
CRP	mg/dl		0.64		0.67		± 0.40	2.44		2.66		± 0.80	4.45		4.85		± 1.50						
	mg/l		6.4		6.7		± 4.0	24.4		26.6		± 8.0	44.5		48.5		± 15.0						

Ref: TEMP-1070 Rev.37 BACK/VERSO 9930089-A