

ABX Minotrol CRP



LOT

MC 401

CONTROL



(Exp.) 2017-01-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio																			
				CONTROL				1		CONTROL				2		CONTROL				3		TOLERANCES TOLERANCE	
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA					
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1		2.1		± 0.4	7.7		8.0		± 1.0	20.5		21.5		± 2.2						
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.43		2.43		± 0.15	4.62		4.61		± 0.18	5.64		5.67		± 0.25						
HB	HGB	g/dl	6.3		6.1		± 0.4	13.6		13.6		± 0.5	17.7		18.1		± 0.7						
		g/l	63		61		± 4	136		136		± 5	177		181		± 7						
		mmol/l	3.91		3.79		± 0.25	8.45		8.45		± 0.31	10.99		11.24		± 0.44						
HT	HCT	%	17.0		17.0		± 2.0	37.4		36.9		± 2.5	49.6		48.8		± 3.0						
		l/l	0.170		0.170		± 0.020	0.374		0.369		± 0.025	0.496		0.488		± 0.030						
VGM	MCV	µm ³ ; fl	70		70.0		± 4	81		80.0		± 4	88		86.0		± 4						
TGMH	MCH	pg	25.9		25.1		± 2.4	29.4		29.5		± 2.4	31.4		31.9		± 2.8						
		fmol	1.61		1.56		± 0.15	1.83		1.83		± 0.15	1.95		1.98		± 0.18						
CCMH	MCHC	g/dl	37.0		35.9		± 3.0	36.3		36.9		± 3.0	35.7		37.1		± 3.0						
		g/l	370		359		± 30	363		369		± 30	357		371		± 30						
		mmol/l	23.00		22.27		± 1.86	22.57		22.90		± 1.86	22.15		23.05		± 1.86						
IDR	RDW	%	13.5		13.5		± 3.0	13.0		13.8		± 3.0	12.5		12.5		± 3.0						
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	72		72		± 20	250		250		± 40	475		498		± 60						
VPM	MPV	µm ³ ; fl	7.5		8.6		± 2.0	7.8		8.6		± 2.0	7.5		7.9		± 2.0						
LYMPHO		#	1.3		N/A		± 0.3	2.6		N/A		± 0.5	3.7		N/A		± 1.3						
		%	63.0		N/A		± 8.0	34.0		N/A		± 6.0	18.0		N/A		± 6.0						
MONO		#	0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0						
		%	9.5		N/A		± 7.0	6.0		N/A		± 5.0	5.0		N/A		± 5.0						
GRANULO		#	0.6		N/A		± 0.3	4.6		N/A		± 0.5	15.8		N/A		± 1.5						
		%	27.5		N/A		± 8.0	60.0		N/A		± 6.0	77.0		N/A		± 7.0						
CRP		mg/dl	0.55		0.53		± 0.40	2.35		2.30		± 0.80	4.05		4.14		± 1.50						
		mg/l	5.5		5.3		± 4.0	23.5		23.0		± 8.0	40.5		41.4		± 15.0						

Ref: TEMP-1070 Rev.36 BACK/VERSO 9930089-A