

# ABX Minotrol CRP



**LOT** MC 072  
Rev 1

**CONTROL**

(Exp.) 2012-10-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio													
				CONTROL			1	CONTROL			2	CONTROL			3		
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE		
				EMICRP				EMICRP				EMICRP					
GB	WBC	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>1.9</b>				± 0.4	<b>7.3</b>				± 1.0	<b>19.4</b>				± 2.2
GR	RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	<b>2.41</b>				± 0.15	<b>4.63</b>				± 0.18	<b>5.51</b>				± 0.25
HB	HGB	g/dl	<b>6.2</b>				± 0.4	<b>13.5</b>				± 0.5	<b>17.7</b>				± 0.7
		g/l	<b>62</b>				± 4	<b>135</b>				± 5	<b>177</b>				± 7
		mmol/l	<b>3.85</b>				± 0.25	<b>8.38</b>				± 0.31	<b>10.99</b>				± 0.44
HT	HCT	%	<b>16.9</b>				± 2.0	<b>38.4</b>				± 2.5	<b>49.6</b>				± 3.0
		l/l	<b>0.169</b>				± 0.020	<b>0.384</b>				± 0.025	<b>0.496</b>				± 0.030
VGM	MCV	µm <sup>3</sup> ; fl	<b>70</b>				± 4	<b>83</b>				± 4	<b>90</b>				± 4
TGMH	MCH	pg	<b>25.7</b>				± 2.4	<b>29.2</b>				± 2.4	<b>32.1</b>				± 2.8
		fmol	<b>1.60</b>				± 0.15	<b>1.81</b>				± 0.15	<b>1.99</b>				± 0.18
CCMH	MCHC	g/dl	<b>36.8</b>				± 3.0	<b>35.1</b>				± 3.0	<b>35.7</b>				± 3.0
		g/l	<b>368</b>				± 30	<b>351</b>				± 30	<b>357</b>				± 30
		mmol/l	<b>22.82</b>				± 1.86	<b>21.82</b>				± 1.86	<b>22.17</b>				± 1.86
IDR	RDW	%	<b>13.0</b>				± 3.0	<b>13.2</b>				± 3.0	<b>12.2</b>				± 3.0
PLAQ.	PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>68</b>				± 20	<b>227</b>				± 40	<b>452</b>				± 60
VPM	MPV	µm <sup>3</sup> ; fl	<b>8.0</b>				± 2.0	<b>7.7</b>				± 2.0	<b>7.6</b>				± 2.0
LYMPHO		#	<b>1.2</b>				± 0.2	<b>2.4</b>				± 0.5	<b>3.0</b>				± 1.0
		%	<b>63.3</b>				± 8.0	<b>33.2</b>				± 6.0	<b>15.7</b>				± 5.0
MONO		#	<b>0.2</b>				± 0.2	<b>0.5</b>				± 0.4	<b>1.0</b>				± 1.0
		%	<b>9.6</b>				± 7.0	<b>6.9</b>				± 5.0	<b>5.1</b>				± 5.0
GRANULO		#	<b>0.5</b>				± 0.2	<b>4.4</b>				± 0.5	<b>15.4</b>				± 1.4
		%	<b>27.0</b>				± 7.0	<b>59.9</b>				± 6.0	<b>79.2</b>				± 7.0
CRP		mg/dl	<b>0.64</b>				± 0.30	<b>2.51</b>				± 0.70	<b>4.21</b>				± 1.00
		mg/l	<b>6.4</b>				± 3.0	<b>25.1</b>				± 7.0	<b>42.1</b>				± 10.0

9930064-A  
Ref: TEMP-1070 Rev.15