

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 043

CONTROL



(Exp.) 2013-07-05

(YYYY - MM - DD)

Rev 2

		ABX Lysebio															
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	
		MICROS					MICROS					MICROS					
		EMI CRP					EMI CRP					EMI CRP					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	1.9				± 0.4	7.3					± 1.0	19.4				± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.47				± 0.15	4.66					± 0.18	5.68				± 0.25
HB HGB	g/dl	5.9				± 0.4	13.3					± 0.5	17.6				± 0.7
	g/l	59				± 4	133					± 5	176				± 7
	mmol/l	3.66				± 0.25	8.26					± 0.31	10.93				± 0.44
HT HCT	%	17.3				± 2.0	39.6					± 2.5	51.1				± 3.0
	l/l	0.173				± 0.020	0.396					± 0.025	0.511				± 0.030
VGM MCV	µm ³ ; fl	70				± 4	85					± 4	90				± 4
TGMH MCH	pg	23.9				± 2.4	28.5					± 2.4	31.0				± 2.8
	fmol	1.48				± 0.15	1.77					± 0.15	1.92				± 0.18
CCMH MCHC	g/dl	34.1				± 3.0	33.6					± 3.0	34.4				± 3.0
	g/l	341				± 30	336					± 30	344				± 30
	mmol/l	21.19				± 1.86	20.85					± 1.86	21.38				± 1.86
IDR RDW	%	13.8				± 3.0	13.8					± 3.0	12.3				± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	73				± 20	261					± 40	504				± 60
VPM MPV	µm ³ ; fl	8.4				± 2.0	8.2					± 2.0	8.1				± 2.0
LYMPHO	#	1.1				± 0.3	2.2					± 0.5	2.8				± 1.1
	%	60.7				± 7.0	30.3					± 6.0	14.2				± 5.0
MONO	#	0.2				± 0.2	0.7					± 0.4	0.9				± 1.0
	%	11.3				± 7.0	9.3					± 5.0	4.9				± 5.0
GRANULO	#	0.6				± 0.3	4.4					± 0.5	15.7				± 1.5
	%	28.1				± 7.0	60.4					± 6.0	80.9				± 7.0
CRP	mg/dl	0.62				± 0.30	2.25					± 0.70	3.90				± 1.00
	mg/l	6.2				± 3.0	22.5					± 7.0	39.0				± 10.0

9930074-A
BACKVERSO
Ref: TEMP-1070 Rev.31