

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 073

CONTROL



(Exp.) 2013-10-05

(YYYY - MM - DD)

Rev 1

		ABX Lysebio															
PARAMETRES PARAMETERS	UNITES UNITS	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	CONTROL				TOLERANCES TOLERANCE	
		MICROS					MICROS					MICROS					
		EMI CRP					EMI CRP					EMI CRP					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0				± 0.4	7.6					± 1.0	19.8				± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.42				± 0.15	4.58					± 0.18	5.52				± 0.25
HB HGB	g/dl	6.1				± 0.4	13.1					± 0.5	17.7				± 0.7
	g/l	61				± 4	131					± 5	177				± 7
	mmol/l	3.79				± 0.25	8.14					± 0.31	10.99				± 0.44
HT HCT	%	16.7				± 2.0	36.6					± 2.5	49.1				± 3.0
	l/l	0.167				± 0.020	0.366					± 0.025	0.491				± 0.030
VGM MCV	µm ³ ; fl	69				± 4	80					± 4	89				± 4
TGMH MCH	pg	25.2				± 2.4	28.6					± 2.4	32.1				± 2.8
	fmol	1.57				± 0.15	1.78					± 0.15	1.99				± 0.18
CCMH MCHC	g/dl	36.5				± 3.0	35.8					± 3.0	36.0				± 3.0
	g/l	365				± 30	358					± 30	360				± 30
	mmol/l	22.69				± 1.86	22.20					± 1.86	22.37				± 1.86
IDR RDW	%	13.1				± 3.0	12.6					± 3.0	11.9				± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	80				± 20	255					± 40	505				± 60
VPM MPV	µm ³ ; fl	8.7				± 2.0	8.2					± 2.0	8.5				± 2.0
LYMPHO	#	1.3				± 0.3	2.7					± 0.5	3.2				± 1.2
	%	63.0				± 8.0	35.0					± 6.0	16.0				± 6.0
MONO	#	0.2				± 0.2	0.5					± 0.4	1.0				± 1.0
	%	9.0				± 7.0	6.5					± 5.0	5.0				± 5.0
GRANULO	#	0.6				± 0.2	4.4					± 0.6	15.6				± 1.4
	%	28.0				± 7.0	58.5					± 6.0	79.0				± 7.0
CRP	mg/dl	0.75				± 0.30	2.47					± 0.70	4.66				± 1.00
	mg/l	7.5				± 3.0	24.7					± 7.0	46.6				± 10.0

9930074-A
BACKVERSO
Ref: TEMP-1070 Rev.31