

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 014

CONTROL



(Exp.) 2014-04-05

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		ABX Lysebio																
		CONTROL				1	CONTROL				2	CONTROL				3		
		MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE		
		EMI CRP					EMI CRP					EMI CRP						
UNITES UNITS																		
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1				± 0.4	7.8					± 1.0	20.4				± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.37				± 0.15	4.60					± 0.18	5.62				± 0.25
HB	HGB	g/dl	6.0				± 0.4	13.3					± 0.5	17.7				± 0.7
		g/l	60				± 4	133					± 5	177				± 7
		mmol/l	3.73				± 0.25	8.26					± 0.31	10.99				± 0.44
HT	HCT	%	16.4				± 2.0	37.3					± 2.5	49.5				± 3.0
		l/l	0.164				± 0.020	0.373					± 0.025	0.495				± 0.030
VGM	MCV	µm ³ , fl	69				± 4	81					± 4	88				± 4
TGMH	MCH	pg	25.3				± 2.4	28.9					± 2.4	31.5				± 2.8
		fmol	1.57				± 0.15	1.80					± 0.15	1.96				± 0.18
CCMH	MCHC	g/dl	36.7				± 3.0	35.7					± 3.0	35.8				± 3.0
		g/l	367				± 30	357					± 30	358				± 30
		mmol/l	22.78				± 1.86	22.17					± 1.86	22.23				± 1.86
IDR	RDW	%	13.2				± 4.0	12.8					± 4.0	11.9				± 4.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	75				± 20	235					± 40	505				± 60
VPM	MPV	µm ³ , fl	8.3				± 2.0	8.5					± 2.0	8.2				± 2.0
LYMPHO		#	1.3				± 0.2	2.6					± 0.5	3.4				± 1.3
		%	63.0				± 8.0	33.5					± 6.0	16.5				± 6.0
MONO		#	0.2				± 0.2	0.5					± 0.4	1.1				± 1.0
		%	9.0				± 7.0	6.0					± 5.0	5.5				± 5.0
GRANULO		#	0.6				± 0.2	4.7					± 0.5	15.9				± 1.5
		%	28.0				± 7.0	60.5					± 6.0	78.0				± 7.0
CRP		mg/dl	0.51				± 0.30	2.28					± 0.70	4.48				± 1.00
		mg/l	5.1				± 3.0	22.8					± 7.0	44.8				± 10.0

Ref: TEMP-1070 Rev.33 BACKVERSO 9930074-A