

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 044

CONTROL



(Exp.) 2014-07-05

(YYYY - MM - DD)

Rev 2

PARAMETRES PARAMETERS		ABX Lysebio																
		CONTROL				1	CONTROL				2	CONTROL				3		
		MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE	MICROS				TOLERANCES TOLERANCE		
		EMI CRP					EMI CRP					EMI CRP						
UNITES UNITS																		
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1				± 0.4	7.7					± 1.0	20.5				± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.43				± 0.15	4.60					± 0.18	5.62				± 0.25
HB	HGB	g/dl	6.1				± 0.4	13.4					± 0.5	17.7				± 0.7
		g/l	61				± 4	134					± 5	177				± 7
		mmol/l	3.79				± 0.25	8.32					± 0.31	10.99				± 0.44
HT	HCT	%	16.8				± 2.0	36.8					± 2.5	50.0				± 3.0
		l/l	0.168				± 0.020	0.368					± 0.025	0.500				± 0.030
VGM	MCV	µm ³ , fl	69				± 4	80					± 4	89				± 4
TGMH	MCH	pg	25.1				± 2.4	29.1					± 2.4	31.5				± 2.8
		fmol	1.56				± 0.15	1.81					± 0.15	1.96				± 0.18
CCMH	MCHC	g/dl	36.4				± 3.0	36.4					± 3.0	35.4				± 3.0
		g/l	364				± 30	364					± 30	354				± 30
		mmol/l	22.59				± 1.86	22.61					± 1.86	21.98				± 1.86
IDR	RDW	%	13.2				± 3.0	12.6					± 3.0	12.1				± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	82				± 20	250					± 40	500				± 60
VPM	MPV	µm ³ , fl	8.6				± 2.0	8.4					± 2.0	8.4				± 2.0
LYMPHO		#	1.4				± 0.3	2.7					± 0.5	3.4				± 1.3
		%	64.5				± 8.0	34.5					± 6.0	16.5				± 6.0
MONO		#	0.2				± 0.2	0.5					± 0.4	1.1				± 1.0
		%	8.5				± 7.0	7.0					± 5.0	5.5				± 5.0
GRANULO		#	0.6				± 0.2	4.5					± 0.5	16.0				± 1.5
		%	27.0				± 7.0	58.5					± 6.0	78.0				± 7.0
CRP		mg/dl	0.50				± 0.30	2.50					± 0.70	4.00				± 1.00
		mg/l	5.0				± 3.0	25.0					± 7.0	40.0				± 10.0

9930074-A BACKVERSO Ref: TEMP-1070 Rev.34