

ABX Minotrol CRP



LOT

MC 400

CONTROL



(Exp.) 2016-10-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio												TOLERANCES TOLERANCE													
				CONTROL				1				CONTROL						2				CONTROL				3			
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA				MICROS		PENTRA		MICROS		PENTRA					
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP				EMI CRP		MS CRP		EMI CRP		MS CRP					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0		2.0		± 0.4	7.7		7.9		± 1.0	20.4		21.2		± 2.2												
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.38		2.37		± 0.15	4.61		4.63		± 0.18	5.70		5.76		± 0.25												
HB	HGB	g/dl	6.2		6.0		± 0.4	13.5		13.5		± 0.5	17.9		18.2		± 0.7												
		g/l	62		60		± 4	135		135		± 5	179		182		± 7												
		mmol/l	3.85		3.73		± 0.25	8.38		8.38		± 0.31	11.12		11.30		± 0.44												
HT	HCT	%	16.7		16.5		± 2.0	37.3		36.3		± 2.5	49.6		49.0		± 3.0												
		l/l	0.167		0.165		± 0.020	0.373		0.363		± 0.025	0.496		0.490		± 0.030												
VGM	MCV	µm ³ ; fl	70		69.5		± 4	81		78.5		± 4	87		85.0		± 4												
TGMH	MCH	pg	26.1		25.3		± 2.4	29.3		29.2		± 2.4	31.4		31.6		± 2.8												
		fmol	1.62		1.57		± 0.15	1.82		1.81		± 0.15	1.95		1.96		± 0.18												
CCMH	MCHC	g/dl	37.2		36.4		± 3.0	36.2		37.1		± 3.0	36.1		37.2		± 3.0												
		g/l	372		364		± 30	362		371		± 30	361		372		± 30												
		mmol/l	23.11		22.62		± 1.86	22.45		23.07		± 1.86	22.42		23.08		± 1.86												
IDR	RDW	%	13.3		13.3		± 3.0	12.4		12.5		± 3.0	12.0		12.6		± 3.0												
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	75		67		± 20	250		245		± 40	500		510		± 60												
VPM	MPV	µm ³ ; fl	8.0		9.1		± 2.0	8.3		8.7		± 2.0	8.3		8.7		± 2.0												
LYMPHO	#		1.3		N/A		± 0.3	2.7		N/A		± 0.5	3.5		N/A		± 1.2												
	%		64.0		N/A		± 8.0	35.5		N/A		± 6.0	17.2		N/A		± 6.0												
MONO	#		0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.1		N/A		± 1.0												
	%		8.0		N/A		± 7.0	6.5		N/A		± 5.0	5.3		N/A		± 5.0												
GRANULO	#		0.6		N/A		± 0.3	4.5		N/A		± 0.5	15.8		N/A		± 1.4												
	%		28.0		N/A		± 8.0	58.0		N/A		± 6.0	77.5		N/A		± 7.0												
CRP	mg/dl		0.47		0.45		± 0.40	2.53		2.58		± 0.80	4.37		4.65		± 1.50												
	mg/l		4.7		4.5		± 4.0	25.3		25.8		± 8.0	43.7		46.5		± 15.0												

Ref: TEMP-1070 Rev.36 BACK/VERSO 9930089-A