

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 015

CONTROL



(Exp.) 2015-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio														
			CONTROL					CONTROL					CONTROL				
			MICROS		1			MICROS		2			MICROS		3		
			EMI CRP		MS CRP		TOLERANCES TOLERANCE	EMI CRP		MS CRP		TOLERANCES TOLERANCE	EMI CRP		MS CRP		TOLERANCES TOLERANCE
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.2		2.1		± 0.4	7.7		7.8		± 1.0	20.6		20.8		± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.39		2.32		± 0.15	4.56		4.52		± 0.18	5.57		5.53		± 0.25
HB	HGB	g/dl	6.2		5.9		± 0.4	13.4		13.5		± 0.5	17.6		18.0		± 0.7
		g/l	62		59		± 4	134		135		± 5	176		180		± 7
		mmol/l	3.85		3.66		± 0.25	8.32		8.38		± 0.31	10.93		11.18		± 0.44
HT	HCT	%	16.5		16.7		± 2.0	36.9		37.1		± 2.5	49.0		48.7		± 3.0
		l/l	0.165		0.167		± 0.020	0.369		0.371		± 0.025	0.490		0.487		± 0.030
VGM	MCV	µm ³ ; fl	69		72.0		± 4	81		82.0		± 4	88		88.1		± 4
TGMH	MCH	pg	25.9		25.4		± 2.4	29.4		29.9		± 2.4	31.6		32.5		± 2.8
		fmol	1.61		1.58		± 0.15	1.82		1.85		± 0.15	1.96		2.02		± 0.18
CCMH	MCHC	g/dl	37.6		35.3		± 3.0	36.3		36.4		± 3.0	35.9		36.9		± 3.0
		g/l	376		353		± 30	363		364		± 30	359		369		± 30
		mmol/l	23.35		21.93		± 1.86	22.53		22.62		± 1.86	22.30		22.94		± 1.86
IDR	RDW	%	13.1		12.0		± 3.0	12.7		11.2		± 3.0	12.3		11.0		± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	79		65		± 20	255		250		± 40	505		510		± 60
VPM	MPV	µm ³ ; fl	8.4		8.9		± 2.0	8.2		8.8		± 2.0	8.0		8.5		± 2.0
LYMPHO		#	1.4		N/A		± 0.2	2.6		N/A		± 0.5	3.3		N/A		± 1.3
		%	64.0		N/A		± 8.0	33.5		N/A		± 6.0	16.0		N/A		± 6.0
MONO		#	0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0
		%	8.5		N/A		± 7.0	7.0		N/A		± 5.0	5.0		N/A		± 5.0
GRANULO		#	0.6		N/A		± 0.2	4.6		N/A		± 0.5	16.3		N/A		± 1.5
		%	27.5		N/A		± 7.0	59.5		N/A		± 6.0	79.0		N/A		± 7.0
CRP		mg/dl	0.50		0.50		± 0.30	2.30		2.20		± 0.70	4.20		4.40		± 1.00
		mg/l	5.0		5.0		± 3.0	23.0		22.0		± 7.0	42.0		44.0		± 10.0

9930089-A
BACKVERSO
Ref: TEMP-1070 Rev.36