

ABX Minotrol CRP



LOT

MC 046

CONTROL



(Exp.) 2016-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio												TOLERANCES TOLERANCE							
				CONTROL			1			CONTROL			2					CONTROL			3		
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA				MICROS		PENTRA			
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP				EMI CRP		MS CRP			
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1		2.1		± 0.4	7.6		7.9		± 1.0	20.1		21.3		± 2.2						
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.31		2.33		± 0.15	4.55		4.61		± 0.18	5.49		5.59		± 0.25						
HB	HGB	g/dl	6.2		6.0		± 0.4	13.3		13.4		± 0.5	17.5		17.9		± 0.7						
		g/l	62		60		± 4	133		134		± 5	175		179		± 7						
		mmol/l	3.85		3.73		± 0.25	8.26		8.32		± 0.31	10.87		11.12		± 0.44						
HT	HCT	%	16.4		16.5		± 2.0	36.9		36.9		± 2.5	48.9		48.6		± 3.0						
		l/l	0.164		0.165		± 0.020	0.369		0.369		± 0.025	0.489		0.486		± 0.030						
VGM	MCV	µm ³ ; fl	71		71.0		± 4	81		80.0		± 4	89		87.0		± 4						
TGMH	MCH	pg	26.8		25.8		± 2.4	29.2		29.1		± 2.4	31.9		32.0		± 2.8						
		fmol	1.67		1.60		± 0.15	1.82		1.81		± 0.15	1.98		1.99		± 0.18						
CCMH	MCHC	g/dl	37.8		36.3		± 3.0	36.1		36.3		± 3.0	35.8		36.8		± 3.0						
		g/l	378		363		± 30	361		363		± 30	358		368		± 30						
		mmol/l	23.48		22.52		± 1.86	22.41		22.56		± 1.86	22.24		22.86		± 1.86						
IDR	RDW	%	13.0		12.5		± 3.0	12.7		13.0		± 3.0	12.3		11.6		± 3.0						
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	70		60		± 20	265		257		± 40	510		510		± 60						
VPM	MPV	µm ³ ; fl	8.7		9.5		± 2.0	8.4		9.0		± 2.0	8.3		8.4		± 2.0						
LYMPHO		#	1.3		N/A		± 0.4	2.7		N/A		± 0.5	3.6		N/A		± 1.3						
		%	63.0		N/A		± 8.0	36.0		N/A		± 6.0	18.0		N/A		± 6.0						
MONO		#	0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0						
		%	9.5		N/A		± 7.0	6.5		N/A		± 5.0	5.0		N/A		± 5.0						
GRANULO		#	0.6		N/A		± 0.3	4.4		N/A		± 0.6	15.5		N/A		± 1.5						
		%	27.5		N/A		± 8.0	57.5		N/A		± 6.0	77.0		N/A		± 7.0						
CRP		mg/dl	0.70		0.70		± 0.40	2.70		2.60		± 0.80	4.30		4.50		± 1.50						
		mg/l	7.0		7.0		± 4.0	27.0		26.0		± 8.0	43.0		45.0		± 15.0						

Ref: TEMP-1070 Rev.36 BACK/VERSO 9930089-A