

ABX Minotrol CRP



LOT

MC 420

CONTROL



(Exp.) 2021-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio														TOLERANCES TOLERANCE		
				CONTROL				1	CONTROL				2	CONTROL						3
				MICROS		PENTRA		TOLERANCES TOLERANCE	MICROS		PENTRA		TOLERANCES TOLERANCE	MICROS		PENTRA				TOLERANCES TOLERANCE
				EMI CRP	LC-767G	MS CRP	EMI CRP		LC-767G	MS CRP	EMI CRP	LC-767G		MS CRP						
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.1	± 0.4	8.0	7.9	8.1	± 1.0	20.7	20.7	21.5	± 2.2						
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.40	2.39	2.34	± 0.15	4.59	4.53	4.52	± 0.18	5.51	5.47	5.52	± 0.25						
HB	HGB	g/dl	6.2	6.1	5.9	± 0.4	13.6	13.5	13.5	± 0.5	17.6	17.6	17.9	± 0.7						
		g/l	62	61	59	± 4	136	135	135	± 5	176	176	179	± 7						
		mmol/l	3.85	3.79	3.66	± 0.25	8.45	8.38	8.38	± 0.31	10.93	10.93	11.12	± 0.44						
HT	HCT	%	16.8	16.7	16.4	± 2.0	37.6	37.1	36.8	± 2.5	49.0	49.2	48.9	± 3.0						
		l/l	0.168	0.167	0.164	± 0.020	0.376	0.371	0.368	± 0.025	0.490	0.492	0.489	± 0.030						
VGM	MCV	µm ³ ; fl	70	70	70.0	± 4	82	82	81.5	± 4	89	90	88.5	± 4						
TGMH	MCH	pg	25.8	25.5	25.2	± 2.4	29.6	29.8	29.9	± 2.4	31.9	32.2	32.4	± 2.8						
		fmol	1.60	1.58	1.57	± 0.15	1.84	1.85	1.85	± 0.15	1.98	2.00	2.01	± 0.18						
CCMH	MCHC	g/dl	36.9	36.5	36.0	± 3.0	36.1	36.3	36.6	± 3.0	35.9	35.8	36.6	± 3.0						
		g/l	369	365	360	± 30	361	363	366	± 30	359	358	366	± 30						
		mmol/l	22.92	22.64	22.37	± 1.86	22.44	22.57	22.76	± 1.86	22.29	22.20	22.75	± 1.86						
IDR	RDW	%	14.0	13.0	15.5	± 3.0	12.5	12.5	13.5	± 3.0	12.0	12.0	12.5	± 3.0						
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	73	79	66	± 20	262	260	248	± 40	490	490	495	± 60						
VPM	MPV	µm ³ ; fl	8.5	9.1	10.1	± 2.0	8.8	8.7	9.3	± 2.0	8.5	8.5	9.0	± 2.0						
LYMPHO		#	1.3	1.3	N/A	± 0.2	3.0	2.8	N/A	± 0.5	3.5	3.3	N/A	± 1.1						
		%	67.0	66.0	N/A	± 8.0	37.0	36.0	N/A	± 6.0	17.0	16.0	N/A	± 5.0						
MONO		#	0.2	0.2	N/A	± 0.2	0.4	0.5	N/A	± 0.4	1.3	1.3	N/A	± 1.3						
		%	5.5	6.5	N/A	± 5.5	5.5	6.0	N/A	± 5.5	6.5	6.5	N/A	± 6.5						
GRANULO		#	0.6	0.6	N/A	± 0.2	4.6	4.6	N/A	± 0.5	15.8	16.0	N/A	± 1.5						
		%	27.5	27.5	N/A	± 7.0	57.5	58.0	N/A	± 6.0	76.5	77.5	N/A	± 7.0						
CRP		mg/dl	0.73	0.83	0.63	± 0.40	2.48	2.90	2.21	± 0.80	4.29	4.88	4.03	± 1.50						
		mg/l	7.3	8.3	6.3	± 4.0	24.8	29.0	22.1	± 8.0	42.9	48.8	40.3	± 15.0						

Ref: TEMP-1070 Rev.39 BACKVERSO 1300092170