

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 014

CONTROL



(Exp.) 2014-04-05

Rev 2

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		Minilysebio												TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3		
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE		
				EMI CRP	CRP 200			EMI CRP	CRP 200			EMI CRP	CRP 200				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1	2.1		± 0.4	7.8	7.6		± 1.0	20.4	19.8		± 2.2				
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.37	2.35		± 0.15	4.60	4.64		± 0.18	5.62	5.72		± 0.25				
HB HGB	g/dl	6.0	5.7		± 0.4	13.3	13.2		± 0.5	17.7	17.8		± 0.7				
	g/l	60	57		± 4	133	132		± 5	177	178		± 7				
	mmol/l	3.73	3.54		± 0.25	8.26	8.22		± 0.31	10.99	11.05		± 0.44				
HT HCT	%	16.4	16.2		± 2.0	37.3	37.6		± 2.5	49.5	50.3		± 3.0				
	l/l	0.164	0.162		± 0.020	0.373	0.376		± 0.025	0.495	0.503		± 0.030				
VGM MCV	µm ³ ; fl	69	69		± 4	81	81		± 4	88	88		± 4				
TGMH MCH	pg	25.3	24.3		± 2.4	28.9	28.5		± 2.4	31.5	31.1		± 2.8				
	fmol	1.57	1.51		± 0.15	1.80	1.77		± 0.15	1.96	1.93		± 0.18				
CCMH MCHC	g/dl	36.7	35.2		± 3.0	35.7	35.2		± 3.0	35.8	35.4		± 3.0				
	g/l	367	352		± 30	357	352		± 30	358	354		± 30				
	mmol/l	22.78	21.83		± 1.86	22.17	21.86		± 1.86	22.23	21.96		± 1.86				
IDR RDW	%	13.2	13.5		± 4.0	12.8	13.0		± 4.0	11.9	12.8		± 4.0				
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	75	79		± 20	235	240		± 40	505	520		± 60				
VPM MPV	µm ³ ; fl	8.3	8.3		± 2.0	8.5	8.0		± 2.0	8.2	7.7		± 2.0				
LYMPHO	#	1.3	1.4		± 0.2	2.6	2.6		± 0.5	3.4	3.1		± 1.3				
	%	63.0	67.6		± 8.0	33.5	34.7		± 6.0	16.5	15.5		± 6.0				
MONO	#	0.2	0.2		± 0.2	0.5	0.5		± 0.4	1.1	1.3		± 1.0				
	%	9.0	8.0		± 7.0	6.0	6.7		± 5.0	5.5	6.5		± 5.0				
GRANULO	#	0.6	0.5		± 0.2	4.7	4.5		± 0.5	15.9	15.5		± 1.5				
	%	28.0	24.5		± 7.0	60.5	58.8		± 6.0	78.0	78.2		± 7.0				
CRP	mg/dl	0.51	0.60		± 0.30	2.28	2.37		± 0.70	4.48	4.49		± 1.00				
	mg/l	5.1	6.0		± 3.0	22.8	23.7		± 7.0	44.8	44.9		± 10.0				

9930087-B
Ref: TEMP-1070 Rev.33