

ABX Minotrol CRP



LOT

MC 420

Rev 1

CONTROL



(Exp.) **2021-10-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio															
				CONTROL				1	CONTROL				2	CONTROL				3	TOLERANCES TOLERANCE
				MICROS		PENTRA		TOLERANCES TOLERANCE	MICROS		PENTRA		TOLERANCES TOLERANCE	MICROS		PENTRA			
				EMI CRP	LC-767G	MS CRP			EMI CRP	LC-767G	MS CRP			EMI CRP	LC-767G	MS CRP			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0		2.1		± 0.4	8.0	7.9		8.1		± 1.0	20.7	20.7		21.5		± 2.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.40	2.39		2.34		± 0.15	4.59	4.53		4.52		± 0.18	5.51	5.47		5.52		± 0.25
HB HGB	g/dl	6.2	6.1		5.9		± 0.4	13.6	13.5		13.5		± 0.5	17.6	17.6		17.9		± 0.7
	g/l	62	61		59		± 4	136	135		135		± 5	176	176		179		± 7
	mmol/l	3.85	3.79		3.66		± 0.25	8.45	8.38		8.38		± 0.31	10.93	10.93		11.12		± 0.44
HT HCT	%	16.8	16.7		16.4		± 2.0	37.6	37.1		36.8		± 2.5	49.0	49.2		48.9		± 3.0
	l/l	0.168	0.167		0.164		± 0.020	0.376	0.371		0.368		± 0.025	0.490	0.492		0.489		± 0.030
VGM MCV	µm ³ ; fl	70	70		70.0		± 4	82	82		81.5		± 4	89	90		88.5		± 4
TGMH MCH	pg	25.8	25.5		25.2		± 2.4	29.6	29.8		29.9		± 2.4	31.9	32.2		32.4		± 2.8
	fmol	1.60	1.58		1.57		± 0.15	1.84	1.85		1.85		± 0.15	1.98	2.00		2.01		± 0.18
CCMH MCHC	g/dl	36.9	36.5		36.0		± 3.0	36.1	36.3		36.6		± 3.0	35.9	35.8		36.6		± 3.0
	g/l	369	365		360		± 30	361	363		366		± 30	359	358		366		± 30
	mmol/l	22.92	22.64		22.37		± 1.86	22.44	22.57		22.76		± 1.86	22.29	22.20		22.75		± 1.86
IDR RDW	%	14.0	13.0		15.5		± 3.0	12.5	12.5		13.5		± 3.0	12.0	12.0		12.5		± 3.0
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	73	79		66		± 20	262	260		248		± 40	490	490		495		± 60
VPM MPV	µm ³ ; fl	8.5	9.1		10.1		± 2.0	8.8	8.7		9.3		± 2.0	8.5	8.5		9.0		± 2.0
LYMPHO	#	1.3	1.3		N/A		± 0.2	3.0	2.8		N/A		± 0.5	3.5	3.3		N/A		± 1.1
	%	67.0	66.0		N/A		± 8.0	37.0	36.0		N/A		± 6.0	17.0	16.0		N/A		± 5.0
MONO	#	0.2	0.2		N/A		± 0.2	0.4	0.5		N/A		± 0.4	1.3	1.3		N/A		± 1.3
	%	5.5	6.5		N/A		± 5.5	5.5	6.0		N/A		± 5.5	6.5	6.5		N/A		± 6.5
GRANULO	#	0.6	0.6		N/A		± 0.2	4.6	4.6		N/A		± 0.5	15.8	16.0		N/A		± 1.5
	%	27.5	27.5		N/A		± 7.0	57.5	58.0		N/A		± 6.0	76.5	77.5		N/A		± 7.0
CRP	mg/dl	0.73	0.83		0.63		± 0.40	2.48	2.90		2.21		± 0.80	4.29	4.88		4.03		± 1.50
	mg/l	7.3	8.3		6.3		± 4.0	24.8	29.0		22.1		± 8.0	42.9	48.8		40.3		± 15.0

Ref: TEMP-1070 Rev.39 BACKVERSO 1300092170