

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 074

CONTROL



(Exp.) 2014-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio																	
			CONTROL			1		TOLERANCES TOLERANCE	CONTROL			2		TOLERANCES TOLERANCE	CONTROL			3		TOLERANCES TOLERANCE
			MICROS		PENTRA	EMI CRP	MICROS		PENTRA	EMI CRP	MICROS		PENTRA		EMI CRP	MICROS		PENTRA		
			EMI CRP		MS CRP		EMI CRP				MS CRP	EMI CRP				MS CRP	EMI CRP		MS CRP	
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1		2.0	± 0.4	7.8		7.8		± 1.0	20.2		20.8		± 2.2				
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.38		2.31	± 0.15	4.62		4.53		± 0.18	5.45		5.37		± 0.25				
HB	HGB	g/dl	6.0		5.7	± 0.4	13.4		13.3		± 0.5	17.5		17.7		± 0.7				
		g/l	60		57	± 4	134		133		± 5	175		177		± 7				
		mmol/l	3.73		3.54	± 0.25	8.32		8.26		± 0.31	10.87		10.99		± 0.44				
HT	HCT	%	16.4		16.5	± 2.0	37.0		36.5		± 2.5	48.5		48.1		± 3.0				
		l/l	0.164		0.165	± 0.020	0.370		0.365		± 0.025	0.485		0.481		± 0.030				
VGM	MCV	µm ³ ; fl	69		71.4	± 4	80		80.6		± 4	89		89.6		± 4				
TGMH	MCH	pg	25.2		24.7	± 2.4	29.0		29.4		± 2.4	32.1		33.0		± 2.8				
		fmol	1.57		1.53	± 0.15	1.80		1.82		± 0.15	1.99		2.05		± 0.18				
CCMH	MCHC	g/dl	36.5		34.6	± 3.0	36.3		36.4		± 3.0	36.1		36.8		± 3.0				
		g/l	365		346	± 30	363		364		± 30	361		368		± 30				
		mmol/l	22.69		21.46	± 1.86	22.51		22.62		± 1.86	22.40		22.84		± 1.86				
IDR	RDW	%	13.5		11.4	± 3.0	12.8		9.8		± 3.0	12.0		9.0		± 3.0				
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	80		72	± 20	258		265		± 40	500		510		± 60				
VPM	MPV	µm ³ ; fl	8.1		9.0	± 2.0	7.8		8.4		± 2.0	8.0		8.4		± 2.0				
LYMPHO		#	1.3		N/A	± 0.2	2.7		N/A		± 0.5	3.4		N/A		± 1.1				
		%	63.5		N/A	± 8.0	34.5		N/A		± 6.0	17.0		N/A		± 5.0				
MONO		#	0.2		N/A	± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0				
		%	8.0		N/A	± 7.0	7.0		N/A		± 5.0	5.0		N/A		± 5.0				
GRANULO		#	0.6		N/A	± 0.2	4.6		N/A		± 0.5	15.8		N/A		± 1.5				
		%	28.5		N/A	± 7.0	58.5		N/A		± 6.0	78.0		N/A		± 7.0				
CRP		mg/dl	0.55		0.50	± 0.30	2.50		2.40		± 0.70	4.10		3.90		± 1.00				
		mg/l	5.5		5.0	± 3.0	25.0		24.0		± 7.0	41.0		39.0		± 10.0				

9930089-A
BACKVERSO
Ref: TEMP-1070 Rev.35