

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 104

CONTROL



(Exp.) 2015-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio												TOLERANCES TOLERANCE	
				CONTROL				CONTROL				CONTROL					
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA			
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP			
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	2.1		2.2		± 0.4	7.8		7.9		± 1.0	20.6		21.4		± 2.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.26		2.20		± 0.15	4.48		4.48		± 0.18	5.48		5.47		± 0.25
HB	HGB	g/dl	6.0		5.8		± 0.4	13.3		13.5		± 0.5	17.7		18.1		± 0.7
		g/l	60		58		± 4	133		135		± 5	177		181		± 7
		mmol/l	3.73		3.60		± 0.25	8.26		8.38		± 0.31	10.99		11.24		± 0.44
HT	HCT	%	16.0		15.6		± 2.0	36.7		35.6		± 2.5	49.3		48.1		± 3.0
		l/l	0.160		0.156		± 0.020	0.367		0.356		± 0.025	0.493		0.481		± 0.030
VGM	MCV	µm ³ , fl	71		70.9		± 4	82		79.5		± 4	90		87.9		± 4
TGMH	MCH	pg	26.5		26.4		± 2.4	29.7		30.1		± 2.4	32.3		33.1		± 2.8
		fmol	1.65		1.64		± 0.15	1.84		1.87		± 0.15	2.01		2.05		± 0.18
CCMH	MCHC	g/dl	37.4		37.2		± 3.0	36.2		37.9		± 3.0	35.9		37.6		± 3.0
		g/l	374		372		± 30	362		379		± 30	359		376		± 30
		mmol/l	23.22		23.09		± 1.86	22.48		23.54		± 1.86	22.29		23.38		± 1.86
IDR	RDW	%	13.0		12.5		± 3.0	12.5		11.4		± 3.0	12.0		11.2		± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	77		65		± 20	255		255		± 40	495		513		± 60
VPM	MPV	µm ³ , fl	8.4		9.4		± 2.0	8.0		9.0		± 2.0	8.1		8.8		± 2.0
LYMPHO		#	1.3		N/A		± 0.2	2.7		N/A		± 0.5	3.4		N/A		± 1.3
		%	63.0		N/A		± 8.0	35.0		N/A		± 6.0	16.5		N/A		± 6.0
MONO		#	0.2		N/A		± 0.2	0.5		N/A		± 0.4	1.0		N/A		± 1.0
		%	9.0		N/A		± 7.0	7.0		N/A		± 5.0	5.0		N/A		± 5.0
GRANULO		#	0.6		N/A		± 0.2	4.5		N/A		± 0.5	16.2		N/A		± 1.5
		%	28.0		N/A		± 7.0	58.0		N/A		± 6.0	78.5		N/A		± 7.0
CRP		mg/dl	0.55		0.50		± 0.30	2.30		2.20		± 0.70	3.90		3.60		± 1.00
		mg/l	5.5		5.0		± 3.0	23.0		22.0		± 7.0	39.0		36.0		± 10.0

9930089-A
BACKVERSO
Ref: TEMP-1070 Rev.36