

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 015

CONTROL



(Exp.) 2015-04-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio												TOLERANCES TOLERANCE													
				CONTROL				1				CONTROL						2				CONTROL				3			
				MICROS		PENTRA		MICROS		PENTRA		MICROS		PENTRA				MICROS		PENTRA		MICROS		PENTRA					
				EMI CRP		MS CRP		EMI CRP		MS CRP		EMI CRP		MS CRP				EMI CRP		MS CRP		EMI CRP		MS CRP					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l		2.2		2.1		± 0.4		7.7		7.8		± 1.0		20.6		20.8		± 2.2									
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l		2.39		2.32		± 0.15		4.56		4.52		± 0.18		5.57		5.53		± 0.25									
HB	HGB	g/dl		6.2		5.9		± 0.4		13.4		13.5		± 0.5		17.6		18.0		± 0.7									
		g/l		62		59		± 4		134		135		± 5		176		180		± 7									
		mmol/l		3.85		3.66		± 0.25		8.32		8.38		± 0.31		10.93		11.18		± 0.44									
HT	HCT	%		16.5		16.7		± 2.0		36.9		37.1		± 2.5		49.0		48.7		± 3.0									
		l/l		0.165		0.167		± 0.020		0.369		0.371		± 0.025		0.490		0.487		± 0.030									
VGM	MCV	µm ³ ; fl		69		72.0		± 4		81		82.0		± 4		88		88.1		± 4									
TGMH	MCH	pg		25.9		25.4		± 2.4		29.4		29.9		± 2.4		31.6		32.5		± 2.8									
		fmol		1.61		1.58		± 0.15		1.82		1.85		± 0.15		1.96		2.02		± 0.18									
CCMH	MCHC	g/dl		37.6		35.3		± 3.0		36.3		36.4		± 3.0		35.9		36.9		± 3.0									
		g/l		376		353		± 30		363		364		± 30		359		369		± 30									
		mmol/l		23.35		21.93		± 1.86		22.53		22.62		± 1.86		22.30		22.94		± 1.86									
IDR	RDW	%		13.1		12.0		± 3.0		12.7		11.2		± 3.0		12.3		11.0		± 3.0									
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l		79		65		± 20		255		250		± 40		505		510		± 60									
VPM	MPV	µm ³ ; fl		8.4		8.9		± 2.0		8.2		8.8		± 2.0		8.0		8.5		± 2.0									
LYMPHO		#		1.4		N/A		± 0.2		2.6		N/A		± 0.5		3.3		N/A		± 1.3									
		%		64.0		N/A		± 8.0		33.5		N/A		± 6.0		16.0		N/A		± 6.0									
MONO		#		0.2		N/A		± 0.2		0.5		N/A		± 0.4		1.0		N/A		± 1.0									
		%		8.5		N/A		± 7.0		7.0		N/A		± 5.0		5.0		N/A		± 5.0									
GRANULO		#		0.6		N/A		± 0.2		4.6		N/A		± 0.5		16.3		N/A		± 1.5									
		%		27.5		N/A		± 7.0		59.5		N/A		± 6.0		79.0		N/A		± 7.0									
CRP		mg/dl		0.50		0.50		± 0.30		2.30		2.20		± 0.70		4.20		4.40		± 1.00									
		mg/l		5.0		5.0		± 3.0		23.0		22.0		± 7.0		42.0		44.0		± 10.0									

9930089-A
BACKVERSO
Ref: TEMP-1070 Rev.36