

ABX Minotrol Retic



LOT RX 427
Rev 1

CONTROL

(Exp.) 2021-03-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|-----------|--------|----------|-------------------------|-----------|--------|----------|-------------------------|-----------|-------|----------|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | | | |
| | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | | 120 RETIC | DX120 | DX Nexus | |
| Méthode Manuelle / Manual RETIC | % | 1.5 | 1.5 | 1.5 | ± 0.9 | 6.0 | 6.0 | 6.0 | ± 2.7 | 11.4 | 11.4 | 11.4 | ± 3.4 |
| RETIC | % | 1.66 | 1.66 | 1.66 | ± 0.83 | 5.06 | 5.06 | 5.06 | ± 1.27 | 10.30 | 10.30 | 10.30 | ± 2.10 |
| RETIC | 10 ⁶ /µl | 0.057 | 0.057 | 0.057 | ± 0.029 | 0.172 | 0.172 | 0.172 | ± 0.043 | 0.348 | 0.348 | 0.348 | ± 0.070 |
| MFI MRI * | % | 32.6 | 32.6 | 32.6 | ± 8.0 | 45.8 | 45.8 | 45.8 | ± 10.0 | 50.5 | 50.5 | 50.5 | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.41 | 3.41 | 3.41 | ± 0.30 | 3.40 | 3.40 | 3.40 | ± 0.30 | 3.38 | 3.38 | 3.38 | ± 0.30 |
| PIC RANGE | - | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 | 27.0 | 27.0 | 27.0 | ± 5.0 |

* For research use only in Canada and USA

Ref: TEMP-0382 Rev.38 FRONT / RECTO 1219930083-A

ABX Minotrol Retic



LOT RX 427
Rev 1

CONTROL

(Exp.) 2021-03-05
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | ABX | | | | | | | | | | | |
|---------------------------------------|---------------------|--------------|--------|--------|-------------------------|--------------|--------|--------|-------------------------|--------------|--|--|-------------------------|
| | | CONTROL 1 | | | TOLERANCES TOLERANCE | CONTROL 2 | | | TOLERANCES TOLERANCE | CONTROL 3 | | | TOLERANCES TOLERANCE |
| | | PENTRA | PENTRA | PENTRA | | PENTRA | PENTRA | PENTRA | | PENTRA | | | |
| Méthode Manuelle / Manual RETIC | % | 1.5 | | | ± 0.9 | 6.0 | | | ± 2.7 | 11.4 | | | ± 3.4 |
| RETIC | % | 1.83 | | | ± 0.92 | 5.16 | | | ± 1.29 | 10.36 | | | ± 2.08 |
| RETIC | 10 ⁶ /µl | 0.062 | | | ± 0.031 | 0.174 | | | ± 0.044 | 0.347 | | | ± 0.070 |
| MFI MRI * | % | 29.3 | | | ± 8.0 | 43.5 | | | ± 10.0 | 47.8 | | | ± 12.0 |
| GR RBC | 10 ⁶ /µl | 3.39 | | | ± 0.30 | 3.37 | | | ± 0.30 | 3.35 | | | ± 0.30 |
| PIC RANGE | - | 8 | | | ± 4 | 8 | | | ± 4 | 8 | | | ± 4 |

Ref: TEMP-0382 Rev.38 BACK / VERSO 1219930083-A

* For research use only in Canada and USA