

ABX Minotrol Retic



LOT

RX 051

CONTROL



(Exp.) 2011-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL		1	TOLERANCES TOLERANCE	CONTROL		2	TOLERANCES TOLERANCE	CONTROL		3	TOLERANCES TOLERANCE
		PENTRA	PENTRA			PENTRA	PENTRA			PENTRA	PENTRA		
		120 RETIC	DX120		120 RETIC	DX120		120 RETIC	DX120				
Méthode Manuelle / Manual RETIC	%	1.3	1.3		± 0.7	6.1	6.1		± 2.4	11.0	11.0		± 3.3
RETIC	%	1.49	1.75		± 0.88	5.65	5.94		± 1.49	10.90	10.90		± 2.20
RETIC	10 ⁹ /µl	0.052	0.060		± 0.030	0.194	0.202		± 0.051	0.374	0.372		± 0.074
MFI MRI *		33.0	33.2		± 8.0	42.8	43.2		± 10.0	47.3	48.3		± 10.0
GR RBC	10 ⁹ /µl	3.47	3.44		± 0.30	3.44	3.40		± 0.30	3.43	3.41		± 0.30
PIC RANGE		25.5	25.5		± 4.5	25.5	25.5		± 4.5	25.5	25.5		± 4.5

9930034-C

Ref: FORM-0382 Rev.07

* For research use only in Canada and USA